


BATRIK Safe-Soak™ Absorbent Floor Mats are designed to absorb fluids from the operating room floor quickly and efficiently, to ensure that the floor beneath remains slip-free, safe and dry during surgical procedures. A skid-resistant backing is available with the Blue Absorbent Mat. All mats are single use and non-sterile.

- absorbs a variety of fluids
- highly absorbent
- conforms to uneven surfaces
- provides quick and easy clean-up
- latex free, disposable, and non-sterile



Safe-Soak™ Absorbent Floor Mats

Item #	Description	Size	Packaging	Absorbency
SS-20-1000	Safe-Soak™ / White extra absorbent mat	28" x 48" (71cm x 122cm)	20 per box	4.5 litres / mat
SS-20-2000	Safe-Soak™ / Blue absorbent mat with skid-resistant backer	28" x 40" (71cm x 102cm)	20 per box	2.5 litres / mat
SS-20-3000	Safe-Soak™ / Yellow absorbent mat	30" x 52" (76cm x 132cm)	15 per box	3.0 litres / mat

BATRIK Safe-Soak™ SPAC floor suction disk, quickly removes fluid waste from the operating room floor through the application of suction. These suction disks have a highly absorbent base material which allows quick fluid wicking and removal when suction is applied. Safe-Soak™ suction disks are latex free and non sterile.

- Easily maneuverable around operating room floor
- highly absorbent and fast wicking
- fluid is captured from bottom through a highly absorbent material
- latex free, disposable, and non-sterile
- opening permits fluid collection that come into contact with its surface
- highly durable and won't be damaged if stepped on
- the opening allows for easy movement using your foot
- the hose is securely affixed on its side to reduce risk of catching equipment



Safe-Soak™ SPAC

Item #	Description	Size	Packaging
SS-21-1000	Safe-Soak™ SPAC	8.625" / 22cm diameter	12 / box