


LAMINARIA

MedGyn's Laminaria is all-natural cervical disposable and sterile dilator made from a sea-grown plant containing no synthetic materials. Laminaria provides gentle, gradual dilation.

Features & Benefits

MedGyn's disposable and sterile Laminaria is all-natural cervical dilator made from a sea-grown plant containing no synthetic materials. Safely used for over a century, Laminaria provides gentle, gradual dilation. Once inserted into cervical canal, Laminaria absorbs liquid, swelling several times its initial dried diameter, gently dilating the cervix.



Indication for use:

Laminaria should be considered whenever cervical dilation and/or softening of the cervix is desired. Some examples are: cervical stenosis, related to dysmenorrhea, considered a possible cause of infertility, resulting from cauterization or conization, placement and removal of intrauterine devices, induction of labor, radium placement, curettage and drainage of uterine cavity, endometrial biopsy, and suction cannula aspiration.

Standard, approx. length 60mm

021003 Laminaria 3mm Small