

ANESTHESIA CARTS

Our Anesthesia Cart is a flexible storage system for the various anesthesia supplies, medications, and instrumentation required in the operating room. Built to last in stainless steel 304 with a #4 polish, it is easy to clean and maintain.

All anesthesia carts are backed by our **1 year warranty** on structural and all components (excluding casters).

KEY FEATURES

SMART DESIGN

- Comes complete with six full extension drawer assemblies in various depths: 3 x 3" (76 mm), 2 x 6" (152 mm), 1 x 9" (228 mm)
- Smooth surface with a three-sided guard rail
- Over bridge assembly with optional accessory bins and/or shelf provides easily accessible additional storage space

SUPERIOR MANEUVERABILITY

- Caster assemblies with high load bearing capacity and grey non-marking rubber tires provide superior maneuverability
- Full perimeter bumper assembly protects cart body during transport

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS (L x W x H)
OR-3680	ANAESTHESIA CART	27" X 22" X 40" (686 mm x 559 mm x 1016 mm)



OR-3680

Products shown in Anaesthesia Cart not included

ACCESSORIES

ITEM #	DESCRIPTION
OR-3680-DV3	DIVIDER KIT for 3" / 76 mm drawer
OR-3680-DV6	DIVIDER KIT for 6" / 152 mm drawer
OR-3680-GH	GLOVE BOX HOLDER
OR-3680-NC3	NARCOTIC CABINET WITH LOCK inside 3" / 76 mm drawer
OR-3680-NC6	NARCOTIC CABINET WITH LOCK inside 6" / 152 mm drawer
OR-3680-OR	OVERHEAD RAIL ASSEMBLY FOR BINS Complete with 2 lateral supports
OR-3680-PB	POWER BAR, HOSPITAL GRADE Complete with 10 ft / 3 M cord
OR-3680-SH	OVERHEAD SHELF ASSEMBLY 26" L x 10" W / 660 mm L x 254 mm W
OR-3680-SR	SIDE RAIL BRACKET
OR-3680-VS	VELCRO STRAP / BRACKET KIT

