

Your First ™

FOAMING SPRAY

READY-TO-USE ENZYME FOAMING SPRAY

Our Foaming Spray is ideal for instrument pre-cleaning as the thick foam spray enzymatically breaks down all bioburden to allow for effective cleaning and sterilization.

Formulated with **PrAmLi Power³ Enzymes™**, a unique blend of three powerful enzymes along with a surfactant to provide the ultimate in cleaning efficacy



INSTRUCTIONS

1. Shake well
2. Spray immediately after procedure
3. Keep wet until the decontamination process is started

PRODUCT SELECTION

ITEM #	PACKAGING	CASE DIMENSIONS
SS-21-4001	1 x 946 ml (32 oz) Spray Bottle	10" x 4"
SS-21-4012	12 x 946 ml (32 oz) Spray Bottles	18" x 12" x 11"

Important: In addition to our enzymatic formulation, our spray bottle contains as much as 25% more solution when compared to the competition.