

# BATRIK ENDOSCOPY BITESTRONG™ BITE BLOCKS

## ADULT & PEDIATRIC BITE BLOCKS

Our Bitestrong™ Bite Blocks provide a safe passage through the patient's mouth during TEE, protecting the patient's teeth and endoscope.

Our molding process ensures a smooth finish with no sharp edges (flashing) that can often irritate a patient's mouth. All Bite Blocks are individually wrapped in sealed, easy to open PE bags and are latex free, lead free, and DEHP free. Additional straps are available, if required

### SINGLE USE

- Available without dental retention

### REUSABLE

- Designed for repeated usage and will withstand commonly practiced chemical and steam sterilization cycles (134°C / 275°F).
- Available with or without an adjustable retention strap

## PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING	CASE DIMENSIONS
<b>B1</b>	Single-Use <a href="#">Pediatric</a>	16 mm / 48FR, includes strap	100 / Box	10" x 10" x 10"
<b>B1-25</b>	Single-Use <a href="#">Pediatric</a>	16 mm / 48FR, includes strap	25 / Box	10" x 7" x 7"
<b>B1-N</b>	Single-Use <a href="#">Pediatric</a>	16 mm / 48FR, without strap	100 / Box	10" x 10" x 10"
<b>B3</b>	Single-Use	20 mm / 60FR, includes strap	100 / Box	10" x 10" x 10"
<b>B3-N</b>	Single-Use	20 mm / 60FR, without strap	100 / Box	10" x 10" x 10"
<b>B3-SA</b>	Single-Use	20 mm / 60FR, with strap <u>attached</u> to one side	100 / Box	10" x 10" x 10"
<b>B3-NDR</b>	Single-Use	20 mm / 60FR, no dental retention rim	100 / Box	10" x 10" x 10"
<b>B3-NDR-N</b>	Single-Use	Same as B3-N, no dental retention rim, without strap.	100 / Box	10" x 10" x 10"
<b>B4</b>	Single-Use	19 mm / 60FR	100 / Box	10" x 10" x 10"
<b>B4T</b>	Single-Use	19 mm / 60FR, includes strap	100 / Box	10" x 10" x 10"
<b>B4T-N</b>	Single-Use	19 mm / 60FR, without strap	100 / Box	10" x 10" x 10"
<b>B6</b>	Reusable	19 mm / 60FR	100 / Box	10" x 10" x 10"
<b>B6T</b>	Reusable	19 mm / 60FR, includes strap	100 / Box	10" x 10" x 10"
<b>B6T-N</b>	Reusable	19 mm / 60FR, without strap	100 / Box	10" x 10" x 10"
<b>B6-25</b>	Reusable	19 mm / 60FR	25 / Box	10" x 7" x 7"
<b>B6T-25</b>	Reusable	19 mm / 60FR, includes strap	25 / Box	10" x 7" x 7"
<b>B6T-N-25</b>	Reusable	19 mm / 60FR, without strap	25 / Box	10" x 7" x 7"
<b>BS</b>	Single-Use	Bite-Block straps 21.7" x 0.79" / 55 cm x 2 cm	100 / Box	10" x 13" x 1.5"



**BS**



**B1**



**B3**



**B4**



**B6**