

# KERRI BEAR MATERNITY BASSINET™

The Kerri Bear Maternity Bassinet™ is easy to handle and clean, making it the ideal choice for hospitals looking for a versatile, safe, stylish, and long-lasting bassinet.

## PRODUCT COMPONENTS

- Infant carrier
- Integrated changing table
- Spacious lower storage cabinet available in 5 colours: Designer White, Silverpine, Featherstone, Shadow Blue, and Mauve Blush

## KEY FEATURES

### ANTIMICROBIAL COATING

- Manufactured with our specialized lead-free, antibacterial protective powder paint coating\*
- Prevents the spread of dangerous micro-organisms and suppresses the growth of mold and bacteria\*\*
- Durable coating resists chipping and discoloration

### BASSINET BASKET

- Frame and carrier are fully welded and painted with a baked protective anti-bacterial powder coating
- Carrier frame accommodates three basket positions: Neutral, Trendelenburg, and Reverse Trendelenburg
- Basket carrier can be easily adjusted back to reveal a convenient changing table underneath (pull safety switch)
- Safety switch mechanism locks basket carrier in place to prevent movement during transport

### BASE CABINET

- Easy-pull handles and magnetic catches for easy door opening and secure door closing
- Fabricated from easy to disinfect high-quality materials (no wood or laminate components used)
- Large storage space for quick access to supplies and materials with optional shelf available



\* Antibacterial paint specifications available upon request  
\*\* Laboratory documentation available upon request

## MATTRESSES

- Two mattresses available: one for basket and one for changing table
- Mattresses made from fire retardant, antibacterial medium density foam material for infant support and comfort
- Changing table features a three-sided guard rail to ensure mattress is secured/ prevent sliding



## CASTERS

- Heavy duty 3" (76 mm) non-marking grey rubber casters facilitate excellent maneuverability of the bassinet
- 2 x Swivel and 2 x Total Lock








CASTERS

## DIMENSIONS (L x W x H)

<b>OVERALL SIZE</b>	32.5" x 19.5" x 37.75" 826 mm x 495 mm x 959 mm	<b>CHANGING TABLE MATTRESS</b>	31" x 18.5" x 1" 787 mm x 470 mm x 25 mm
<b>ACRYLIC BASKET</b>	27.5" x 13.5" x 9.5" (Head End) 700 mm x 340 mm x 241 mm	<b>BASKET MATTRESS</b>	26" x 12.5" x 2" 660 mm x 318 mm x 50 mm
<b>STORAGE CABINET (INTERIOR)</b>	31" x 18" x 20" 787 mm x 457 mm x 508 mm		

## PRODUCT SELECTION

ITEM #	COLOUR	
PH-2695WH	DESIGNER WHITE	
PH-2695SP	SILVERPINE	
PH-2695FS	FEATHERSTONE	
PH-2695SB	SHADOW BLUE	
PH-2695MB	MAUVE BLUSH	

## ACCESSORIES

ITEM #	DESCRIPTION (L x W x H)
1 PH-2657	MATERNITY BASKET WITH CARD HOLDER High impact seamless acrylic clear Size: 27.5" x 13.5" x 9.5" / 698.5 mm x 342.9 mm x 241.3 mm Card Holder Size: 5" x 3.5
2 PH-2659	BASKET MATTRESS Fire retardant anti-bacterial Size: 26" x 12.5" x 2" / 660.4 mm x 317.5 mm x 50.8 mm
3 PH-2658	CHANGING TABLE MATTRESS Fire retardant anti-bacterial Size: 31" x 18.5" x 1" / 787.4 mm x 469.9 mm x 25.4 mm
4 PH-2672	INTERIOR CABINET WITH OPTIONAL SHELF Not adjustable
5 12008	CASTER 3" ALL SWIVEL All surface
5 12009	CASTER 3" WHEEL LOCK All surface, wheel lock
6 PH-2660	JAR RACK White only
7 PH-2661D	BASIN HOLDER, AND STAINLESS STEEL BASIN
8 PH-2671	CHART HOLDER

# Find the images for the Items codes through the document



FOR EVERY NEW KERRI BEAR  
BASSINET™ SOLD,  
A DONATION OF \$10 WILL  
BE MADE TO THE MISSING  
CHILDREN'S NETWORK<sup>3</sup>.

3. Imperial Surgical has established special arrangements with the respective charity.

## WE OFFER ONE OF THE BEST WARRANTIES IN THE INDUSTRY

All Imperial Surgical Kerri Bear Bassinets™ are warranted to be free from defects in workmanship and materials for a period of **seven years from date of original purchase**.

Accessories which are not manufactured by Imperial Surgical are warranted from defects in workmanship for one year.

