

BATRIK PrAmLi Power³ LUBE™ Instrument Lubricant

Specially formulated to lubricate surgical instruments, hinges, box locks, and general equipment during autoclaving. PrAmLi Power³ LUBE™ is uniquely formulated with an all natural corrosion inhibitor.

- Forms a protective coating which can extend the instrument's, useful life
- Bacteriostatic
- Neutral pH 6.5 - 7.5
- Prevents corrosion

TIP:
FOR MILK BATH SOLUTION,
DILUTE 1 PART SOLUTION
IN 6 PARTS WATER

Formulated with
PrAmLi Power³ Enzymes™

ITEM #	DIMENSIONS
OP-11	0.25" x 0.75" (6.35 mm x 19.05 mm)
OP-12	0.375" x 1" (9.53 mm x 25.4 mm)
OP-13	0.5" x 1" (12.7 mm x 25.4 mm)
OP-14	0.75" x 1" (19.05 mm x 25.4 mm)
OP-15	1" x 1" (25.4 mm x 25.4 mm)
OP-16	1.25" x 1.25" (31.75 mm x 31.75 mm)
OP-17	1.5" x 1.5" (38.1 mm x 38.1 mm)
OP-18	1.75" x 1.25" (44.45 mm x 31.75 mm)
OP-20	Assorted



PACKAGING: Bag of 50

PRODUCT SELECTION

ITEM #	DILUTION LUBRICANT / WATER	PACKAGING	CASE DIMENSIONS
SS-21-6001	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	1 x 3.8 L (1 Gallon)	6" x 6" x 13"
SS-21-6004	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	4 x 3.8 L (4 Gallons)	13" x 13" x 14"
SS-21-6025	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	2 x 9.5 L (2 x 2.5 Gallons)	14.5" x 9.25" x 15"
SS-21-6032	No dilution required	1 x 946 ml (32 oz) Spray bottle	10" x 4"
SS-21-6012	No dilution required	12 X 946 ml (32 oz) Spray Bottles	18" x 12" x 11"