

Your First ™

TRANSPORT BAGS FOR ENDOSCOPE AND MEDICAL INSTRUMENTS

Thoughtfully designed to eliminate touch-points during endoscope transport from point-of-care to the reprocessing area, and back, focused on safety of patients and medical staff.

- Made of recyclable LDPE material known to be quite flexible, highly resilient, and tough
- Latex and DEHP-free
- Excellent resistance to dilute and concentrated acids, alcohols, bases and esters*
- Good resistance to aldehydes, ketones, and vegetable-based oils*
- Can withstand high temperatures of up to 90 °C (194°F)
- Tamper-evident seal for added safety, also aids in keeping the instrument moist to avoid drying of bio-debris
- Clearly marked bags with color-coded visuals facilitates immediate identification of clean vs. used instruments
- Users can write important information directly on the bag in spaces marked with visual icons
- Available in two convenient sizes

GOOD TO KNOW: USING **Your First ™** SEALABLE TRANSPORT BAGS REDUCES USER ERROR, PREVENTS CROSS-CONTAMINATION, AND PROMOTES PATIENT AND MEDICAL STAFF SAFETY.

PRODUCT SELECTION


| ITEM # | DESCRIPTION | DIMENSIONS | PACKAGING | CASE DIMENSIONS |
|-----------------|------------------------|---------------------------------|------------|-----------------|
| 17-HAP24 | Clean Instruments Bag | 24" x 22" / 61.0 cm x 56.0 cm | 200 / Case | 18" x 14" x 5" |
| 17-SAD24 | Soiled Instruments Bag | 24" x 22" / 61.0 cm x 56.0 cm | 200 / Case | 18" x 14" x 5" |
| 17-HAP30 | Clean Instruments Bag | 30.5" x 26" / 77.0 cm x 66.0 cm | 200 / Case | 21" x 19" x 7" |
| 17-SAD30 | Soiled Instruments Bag | 30.5" x 26" / 77.0 cm x 66.0 cm | 200 / Case | 21" x 19" x 7" |

*Reference on file

CLEAN INSTRUMENTS BAG



SOILED INSTRUMENTS BAG



MADE IN CANADA