

BRUSHES PORTFOLIO

Established in 1997. Proudly Canadian. 
Sold in over 85 countries worldwide.

BATRIK ARC™ WIRE STYLE BRUSHES

SINGLE-USE AND REUSABLE

You've got options with our ARC™ Wire Brushes. Our single-use and reusable brushes are available in both flexible and rigid varieties and come in a multitude of sizes. ARC™ Wire Brushes are guaranteed to last longer than twisted wire and preserve the integrity of your valuable instruments.

BRUSH SELECTION TIPS:

INSTRUMENT SPECIFICATIONS

What is the maximum length of the instrument(s) to be cleaned?

Add at least 2" (5cm) to that length, so you have an adequate length to work with.

What is the inner diameter of the instrument?

Add at most 0.125" (3mm) to determine the diameter of the bristle to ensure effective cleaning. The bristle diameter has to slightly exceed the inner tube diameter. Bristles that

are too large will bend too much and not clean effectively.

BRUSH SELECTION

Simply find the section in our catalogue that best corresponds to the overall length and diameter of your instrument. As a rule, if your present brush is 8" (20.32cm) in length and we only offer 10" (25.4cm) length with your required diameter, that is entirely acceptable.

In many cases, our ARC™ Wire Brushes will have more bristles on the brush. As an example, many 5mm (0.19") diameter brushes are only 1" (2.5cm) in length. At BATRIK, we strived to provide the greatest amount of bristles to ensure the most effective cleaning.

CUSTOMIZABLE WIRE BRUSHES


A unique aspect of single-use ARC™ Wire Brushes is that they can be customized to any length:

- Up to 2000mm for Wire Brushes on plastic semi-rigid sheaths, rigid brushes from 5 to 20mm in diameter.
- No limit in length for Wire Brushes on flexible plastic sheaths, flexible brushes from 2-5mm in diameter

RIGID



FLEXIBLE



12" / 30 CM SINGLE-END WIRE BRUSH

| | ITEM # SINGLE-USE | ITEM # REUSABLE | COLOR | TOTAL LENGTH | BRUSH (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|-----------------------------------|----------------------|--------------------|--------|-------------------|-----------------|----------------|-----------------|
| FLEXIBLE BRUSH & SHEATH | 41-31502-S | 41-31502 | White | 12.401" / 31.5 cm | 0.078" / 2 mm | 0.59" / 15 mm | 0.39" / 1 mm |
| | 41-32503-S | 41-32503 | Yellow | 12.795" / 32.5 cm | 0.118" / 3 mm | 0.98" / 25 mm | 0.39" / 1 mm |
| | 411-32504-S | 411-32504 | Violet | 12.795" / 32.5 cm | 0.157" / 4 mm | 0.98" / 25 mm | 0.61" / 1.55 mm |
| | 411-32505-S | 411-32505 | Grey | 12.795" / 32.5 cm | 0.19" / 5 mm | 0.98" / 25 mm | 0.68" / 1.75 mm |
| RIGID BRUSH & SHEATH | 41-35005-S | 41-35005 | White | 13.779" / 35 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| | 41-35006-S | 41-35006 | Green | 13.779" / 35 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| | 41-35008-S | 41-35008 | Blue | 13.779" / 35 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-35010-S | 41-35010 | Red | 13.779" / 35 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-35012-S | 41-35012 | Yellow | 13.779" / 35 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-35015-S | 41-35015 | Green | 13.779" / 35 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| | 41-35020-S | 41-35020 | Blue | 13.779" / 35 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING Single-Use: Bag of 20 / Reusable: Bag of 3

18" / 45 CM SINGLE-END WIRE BRUSH

| | ITEM # SINGLE-USE | ITEM # REUSABLE | COLOR | TOTAL LENGTH | BRUSH (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|-----------------------------------|----------------------|--------------------|--------|-------------------|-----------------|----------------|-----------------|
| FLEXIBLE BRUSH & SHEATH | 41-46502-S | 41-46502 | White | 18.307" / 46.5 cm | 0.078" / 2 mm | 0.59" / 15 mm | 0.39" / 1 mm |
| | 41-47503-S | 41-47503 | Yellow | 18.7" / 47.5 cm | 0.118" / 3 mm | 0.98" / 25 mm | 0.39" / 1 mm |
| | 411-47504-S | 411-47504 | Violet | 18.7" / 47.5 cm | 0.157" / 4 mm | 0.98" / 25 mm | 0.61" / 1.55 mm |
| | 411-47505-S | 411-47505 | Grey | 18.7" / 47.5 cm | 0.19" / 5 mm | 0.98" / 25 mm | 0.68" / 1.75 mm |
| RIGID BRUSH & SHEATH | 41-50005-S | 41-50005 | White | 19.685" / 50 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| | 41-50006-S | 41-50006 | Green | 19.685" / 50 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| | 41-50008-S | 41-50008 | Blue | 19.685" / 50 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-50010-S | 41-50010 | Red | 19.685" / 50 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-50012-S | 41-50012 | Yellow | 19.685" / 50 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-50015-S | 41-50015 | Green | 19.685" / 50 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| | 41-50020-S | 41-50020 | Blue | 19.685" / 50 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Single -Use: Bag of 20 / Reusable: Bag of 3

24" / 60 CM SINGLE-END WIRE BRUSH



| | ITEM # SINGLE-USE | ITEM # REUSABLE | COLOR | TOTAL LENGTH | BRUSH (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|-----------------------------------|----------------------|--------------------|--------|-------------------|-----------------|----------------|-----------------|
| FLEXIBLE BRUSH & SHEATH | 41-61502-S | 41-61502 | White | 24.212" / 61.5 cm | 0.078" / 2 mm | 0.59" / 15 mm | 0.39" / 1mm |
| | 41-62503-S | 41-62503 | Yellow | 24.6" / 62.5 cm | 0.118" / 3 mm | 0.98" / 25 mm | 0.39" / 1mm |
| | 411-62504-S | 411-62504 | Violet | 24.6" / 62.5 cm | 0.157" / 4 mm | 0.98" / 25 mm | 0.61" / 1.55mm |
| | 411-62505-S | 411-62505 | Grey | 24.6" / 62.5 cm | 0.19" / 5 mm | 0.98" / 25 mm | 0.68" / 1.75mm |
| RIGID BRUSH & SHEATH | 41-65005-S | 41-65005 | White | 25.59" / 65 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3mm |
| | 41-65006-S | 41-65006 | Green | 25.59" / 65 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3mm |
| | 41-65008-S | 41-65008 | Blue | 25.59" / 65 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-65010-S | 41-65010 | Red | 25.59" / 65 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-65012-S | 41-65012 | Yellow | 25.59" / 65 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-65015-S | 41-65015 | Green | 25.59" / 65 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| | 41-65020-S | 41-65020 | Blue | 25.59" / 65 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

30" / 75 CM SINGLE-END WIRE BRUSH


| | ITEM # SINGLE-USE | ITEM # REUSABLE | COLOR | TOTAL LENGTH | BRUSH (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|-----------------------------------|----------------------|--------------------|--------|-------------------|-----------------|----------------|-----------------|
| FLEXIBLE BRUSH & SHEATH | 41-76502-S | 41-76502 | White | 30.118" / 76.5 cm | 0.078" / 2 mm | 0.59" / 15 mm | 0.39" / 1 mm |
| | 41-77503-S | 41-77503 | Yellow | 30.511" / 77.5 cm | 0.118" / 3 mm | 0.98" / 25 mm | 0.39" / 1 mm |
| | 411-77504-S | 411-77504 | Violet | 30.511" / 77.5 cm | 0.157" / 4 mm | 0.98" / 25 mm | 0.61" / 1.55 mm |
| | 411-77505-S | 411-77505 | Grey | 30.511" / 77.5 cm | 0.19" / 5 mm | 0.98" / 25 mm | 0.68" / 1.75 mm |
| RIGID BRUSH & SHEATH | 41-80005-S | 41-80005 | White | 31.496" / 80 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| | 41-80006-S | 41-80006 | Green | 31.496" / 80 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| | 41-80008-S | 41-80008 | Blue | 31.496" / 80 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-80010-S | 41-80010 | Red | 31.496" / 80 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-80012-S | 41-80012 | Yellow | 31.496" / 80 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| | 41-80015-S | 41-80015 | Green | 31.496" / 80 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| | 41-80020-S | 41-80020 | Blue | 31.496" / 80 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3



FLEXIBLE

Flexible brush diameter (Ø) ranging from 2 mm to 5 mm



RIGID

Rigid brush diameter (Ø) ranging from 5 mm to 20 mm

BRUSH SELECTION SYSTEM (Actual Size)

| 2mm 0.078 | 3 mm 0.118" | 4 mm 0.157" | 5 mm 0.196" | 6 mm 0.236" | 8 mm 0.314" | 10 mm 0.394" | 12 mm 0.472" | 15 mm 0.591" | 20 mm 0.787" |
|--------------|----------------|----------------|----------------|----------------|----------------|-----------------|-----------------|-----------------|-----------------|
| | | | | | | | | | |

BATRIK REUSABLE DOUBLE ARC™ WIRE BRUSHES

ASYMMETRIC

14" / 35 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH 2 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|-----------------|---------------|-----------------|----------------|-----------------|
| 40-35056 | White | 13.779" / 35 cm | 0.19" / 5 mm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-35810 | Blue | 13.779" / 35 cm | 0.31" / 8 mm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-35112 | Yellow | 13.779" / 35 cm | 0.39" / 10 mm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-35152 | Green | 13.779" / 35 cm | 0.59" / 15 mm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

20" / 50 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH 2 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|-----------------|---------------|-----------------|----------------|-----------------|
| 40-50056 | White | 19.685" / 50 cm | 0.19" / 5 mm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-50810 | Blue | 19.685" / 50 cm | 0.31" / 8 mm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-50112 | Yellow | 19.685" / 50 cm | 0.39" / 10 mm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-50152 | Green | 19.685" / 50 cm | 0.59" / 15 mm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

26" / 65 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH 2 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|----------------|---------------|-----------------|----------------|-----------------|
| 40-65056 | White | 25.59" / 65 cm | 0.19" / 5 mm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-65810 | Blue | 25.59" / 65 cm | 0.31" / 8 mm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-65112 | Yellow | 25.59" / 65 cm | 0.39" / 10 mm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-65152 | Green | 25.59" / 65 cm | 0.59" / 15 mm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

31" / 80 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH 2 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|-----------------|---------------|-----------------|----------------|-----------------|
| 40-80056 | White | 31.496" / 80 cm | 0.19" / 5 mm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-80810 | Blue | 31.496" / 80 cm | 0.31" / 8 mm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-80112 | Yellow | 31.496" / 80 cm | 0.39" / 10 mm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-80152 | Green | 31.496" / 80 cm | 0.59" / 15 mm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

BATRIK REUSABLE DOUBLE ARC™ WIRE BRUSHES SYMMETRIC

14" / 35 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|-----------------|-----------------|----------------|-----------------|
| 40-35005 | White | 13.779" / 35 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-35006 | Green | 13.779" / 35 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-35008 | Blue | 13.779" / 35 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-35010 | Red | 13.779" / 35 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-35012 | Yellow | 13.779" / 35 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-35015 | Green | 13.779" / 35 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| 40-35020 | Blue | 13.779" / 35 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

20" / 50 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|-----------------|-----------------|----------------|-----------------|
| 40-50005 | White | 19.685" / 50 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-50006 | Green | 19.685" / 50 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-50008 | Blue | 19.685" / 50 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-50010 | Red | 19.685" / 50 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-50012 | Yellow | 19.685" / 50 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-50015 | Green | 19.685" / 50 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| 40-50020 | Blue | 19.685" / 50 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

26" / 65 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|----------|--------|----------------|-----------------|----------------|-----------------|
| 40-65005 | White | 25.59" / 65 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-65006 | Green | 25.59" / 65 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-65008 | Blue | 25.59" / 65 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-65010 | Red | 25.59" / 65 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-65012 | Yellow | 25.59" / 65 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-65015 | Green | 25.59" / 65 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| 40-65020 | Blue | 25.59" / 65 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

31" / 80 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

| ITEM # | COLOR | TOTAL LENGTH | BRUSH 1 (Ø) | BRUSH LENGTH | SHEATH (Ø) |
|-----------------|--------|-----------------|-----------------|----------------|-----------------|
| 40-80005 | White | 31.496" / 80 cm | 0.19" / 5 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-80006 | Green | 31.496" / 80 cm | 0.23" / 6 mm | 1.968" / 50 mm | 0.118" / 3 mm |
| 40-80008 | Blue | 31.496" / 80 cm | 0.31" / 8 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-80010 | Red | 31.496" / 80 cm | 0.39" / 10 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-65112 | Yellow | 31.496" / 80 cm | 0.49" / 12.5 mm | 1.968" / 50 mm | 0.157" / 4 mm |
| 40-80015 | Green | 31.496" / 80 cm | 0.59" / 15 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |
| 40-80020 | Blue | 31.496" / 80 cm | 0.78" / 20 mm | 1.968" / 50 mm | 0.228" / 5.8 mm |

PACKAGING: Bag of 3

ASYMMETRIC BRUSHES

5 MM
6 MM

8 MM
10 MM

10 MM
12.5 MM

15 MM
20 MM



BATRIK

REUSABLE IN-BRUSH™ ANTIMICROBIAL TWISTED WIRE BRUSHES

Our In-Brush™ Antimicrobial Reusable Twisted Wire Brushes are designed for use with Frazier/Baron suction tubes. The twisted wires have a fine/medium diameter yet are rigid enough to ensure excellent cleaning capabilities.

Our wide assortment of Nylon Twisted Wire Brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

GOOD TO KNOW: OUR TWISTED WIRE BRUSHES FEATURE A PROTECTIVE BLUE TIP FOR ADDED SAFETY.



10" / 250 MM FRAZIER /BARON SUCTION TUBE BRUSHES

| ITEM # | SIZE | BRUSH LENGTH | BRUSH (Ø) |
|-----------------|-------|----------------|------------------|
| 21-10039 | 3 FR | 0.5" / 12.7 mm | 0.039" / 1 mm |
| 21-10067 | 5 FR | 0.5" / 12.7 mm | 0.065" / 1.67 mm |
| 20-10059 | 5 FR | 1" / 25.4 mm | 0.059" / 1.50 mm |
| 20-10079 | 6 FR | 1" / 25.4 mm | 0.079" / 2 mm |
| 20-10090 | 7 FR | 1" / 25.4 mm | 0.090" / 2.3 mm |
| 20-10106 | 8 FR | 1" / 25.4 mm | 0.106" / 2.7 mm |
| 20-10118 | 9 FR | 1" / 25.4 mm | 0.118" / 3 mm |
| 20-10129 | 10 FR | 1.5" / 38 mm | 0.129" / 3.3 mm |
| 20-10146 | 11 FR | 1.5" / 38 mm | 0.0146" / 3.7 mm |
| 20-10157 | 12 FR | 1.5" / 38 mm | 0.157" / 4 mm |

PACKAGING: Bag of 3.
 Also available in bags of 50, add the suffix "-50".

12" / 305 MM NYLON TWISTED WIRE BRUSHES

| ITEM # | BRUSH LENGTH | BRUSH (Ø) |
|-----------------|--------------|----------------|
| 22-12078 | 2" / 50.8 mm | 0.078" / 2 mm |
| 22-12118 | 2" / 50.8 mm | 0.118" / 3 mm |
| 22-12157 | 3" / 76 mm | 0.157" / 4 mm |
| 22-12197 | 3" / 76 mm | 0.197" / 5 mm |
| 22-12236 | 3" / 76 mm | 0.236" / 6 mm |
| 22-12393 | 3" / 76 mm | 0.394" / 10 mm |

PACKAGING: Bag of 3.
 Also available in bags of 50, add the suffix "- 50".

16" / 406 MM NYLON TWISTED WIRE BRUSHES

| ITEM # | BRUSH LENGTH | BRUSH (Ø) |
|-----------------|----------------|----------------|
| 23-16197 | 2.5" / 63.5 mm | 0.197" / 5 mm |
| 23-16236 | 2.5" / 63.5 mm | 0.236" / 6 mm |
| 23-16276 | 3" / 76 mm | 0.276" / 7 mm |
| 23-16394 | 3" / 76 mm | 0.394" / 10 mm |
| 23-16472 | 3" / 76 mm | 0.472" / 12 mm |
| 23-16591 | 3" / 76 mm | 0.591" / 15 mm |
| 23-16787 | 3" / 76 mm | 0.787" / 20 mm |

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "-50".

18" / 457 MM NYLON TWISTED WIRE BRUSHES

| ITEM # | BRUSH LENGTH | BRUSH (Ø) |
|-----------------|---------------|----------------|
| 19-18010 | 0.47" / 12 mm | 0.039" / 1 mm |
| 19-18015 | 0.5" / 12 mm | 0.06" / 1.6 mm |

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "-50".

24" / 610 MM NYLON TWISTED WIRE BRUSHES


| ITEM # | BRUSH LENGTH | BRUSH (Ø) |
|-----------------|----------------|-----------------|
| 24-24098 | 2.5" / 63.5 mm | 0.098" / 2.5 mm |
| 24-24197 | 2.5" / 63.5 mm | 0.197" / 5 mm |
| 24-24276 | 3" / 76 mm | 0.276" / 7 mm |
| 24-24394 | 3" / 76 mm | 0.394" / 10 mm |
| 24-24472 | 3" / 76 mm | 0.472" / 12 mm |

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "-50".

TRACHEOSTOMY BRUSHES

| ITEM # | DESCRIPTION | TOTAL LENGTH | BRUSH LENGTH | BRUSH (Ø) |
|-----------------|----------------------|--------------|---------------|----------------|
| 25-08500 | Large 6-12.3 3 / 4" | 8" / 203 mm | 3.75" / 95 mm | 0.5" / 12.7 mm |
| 25-37550 | Small 000-5.1 7 / 8" | 4" / 93 mm | 1.88" / 48 mm | 0.5" / 12.7 mm |

PACKAGING: Bag of 3. Also available in bags of 50, add the suffix "-50".



BRONCHO BRUSHES

| ITEM # | DESCRIPTION | BRUSH LENGTH | BRUSH (Ø) INCHES / MM |
|-----------------|-------------|--------------|-----------------------|
| 26-06175 | Small | 6" / 152 mm | 1.75" / 44.5 mm |
| 26-08325 | Large | 8" / 203 mm | 3.25" / 82.6 mm |

PACKAGING: Bag of 3.
Also available in bags of 50, add the suffix "-50".

BATRIK STAINLESS STEEL ARC™ MICRO BRUSHES

Our Stainless Steel ARC™ Micro Brushes are some of the smallest brushes in the industry. They are designed to clean channels with an inner diameter ranging from 0.5 mm to 1.40 mm.

- Flexible stainless steel sheath will not bend like the regular metal twisted wire brushes
- Polyamide bristles prevents or slows bacterial growth between sterilizations
- Reusable*
- Available in two lengths

TIP: USING THE PULL-THROUGH
TECHNIQUE IS HIGHLY
RECOMMENDED WITH THESE
BRUSHES. SIMPLY INSERT THE
BLUNT, NON-BRUSH END FIRST.



PRODUCT SELECTION

| ITEM # | TOTAL LENGTH | BRUSH SIZE (Ø) | BRUSH LENGTH | SHEATH (Ø) | PACKAGING |
|-----------------|--------------|----------------|--------------|-----------------|-----------|
| 35-03015 | 12" / 30 cm | 0.06" / 1.5 mm | 0.6" / 15 mm | 0.02" / 0.40 mm | 6 / Tube |
| 35-04515 | 18" / 45 cm | 0.06" / 1.5 mm | 0.6" / 15 mm | 0.02" / 0.40 mm | 6 / Tube |

*Although all materials used in the fabrication of these brushes can be re-processed repeatedly, due to the delicate nature of these brushes re-use is conditional on their handling and application.

BATRIK INSTRUMENT CLEANING FLAT BRUSHES

Our Instrument Cleaning Flat Brushes are specifically designed for hospital use in the pre-clean and preparation of instruments.

- Long ergonomic handle permits ease-of-use
- Extra long length of brush head (60 [±5] mm in Single-Use and 75 [±5] mm in Reusable)
- Transparent brush fibers allow for quick verification of its cleanliness
- The plastic handle will not release any toxic components
- Available in 6 models (4 reusable/autoclavable and 2 single-use)



FLEXIBLE

This brush is designed for large instruments that require a more gentle cleaning application such as optics, plastics, or hard to reach spot.



EXTRA FLEXIBLE

The very long, flexible bristles eliminate any stress during sanitation, and makes this the ideal brush for fragile instruments as well as both pliable and rigid endoscopes. The length of the bristles assist in cleaning areas that would otherwise be difficult to reach, such as the accessory channel or tip deflection of pliable endoscopes.



RIGID

The rigidity of this brush, along with the length of the bristles and handle, improves and simplifies instrument cleaning.



EXTRA RIGID

Instrument safe polyamide fibers offer steel-like rigidity and allows deep cleaning of soiled instruments, particularly in grooved or uneven surfaces.



50-40919



50-40921



50-40920



50-40918



50-40922



50-40917

PRODUCT SELECTION


| ITEM # | BRUSH TYPE | DESCRIPTION | COLOR | BRUSH LENGTH | BRUSH SIZE | PACKAGING |
|----------|----------------|----------------------------------|------------|--------------|-----------------------|-----------|
| 50-40917 | Extra Rigid | Reusable Autoclavable Flat Brush | Blue | 220 mm | 75 mm x 20 mm x 10 mm | Bag of 3 |
| 50-40918 | Rigid | Reusable Autoclavable Flat Brush | Green | 220 mm | 75 mm x 20 mm x 15 mm | Bag of 3 |
| 50-40919 | Flexible | Reusable Autoclavable Flat Brush | Yellow | 220 mm | 75 mm x 20 mm x 15 mm | Bag of 3 |
| 50-40920 | Extra Flexible | Reusable Autoclavable Flat Brush | Red | 220 mm | 75 mm x 20 mm x 25 mm | Bag of 3 |
| 50-40921 | Flexible | Single-Use Flat Brush | Orange | 220 mm | 60 mm x 20 mm x 15 mm | Bag of 20 |
| 50-40922 | Rigid | Single-Use Flat Brush | Light Blue | 220 mm | 60 mm x 20 mm x 15 mm | Bag of 20 |

Reusable models are also available in bags of 50. Add the suffix "- 50".

BATRIK REUSABLE ORTHOPEDIC RASPS CLEANING BRUSH

Our Reusable Orthopedic Rasps Cleaning Brush is used to clean orthopedic rasps and other similar instruments used in orthopedic surgery.

- Long lasting and can withstand many sterilization cycles
- Ergonomic steamable thermoplastic handle
- Antibacterial polyamide bristles
- Reusable
- Latex free



PRODUCT DESCRIPTION

| ITEM # | TOTAL LENGTH | HANDLE LENGTH | BRUSH SIZE | PACKAGING |
|----------------|--------------|---------------|---|-----------|
| 50-RASP | 9" / 23 cm | 5.5" / 14 cm | 0.78" x 2.7" x 3.5" 2 cm x 7 cm x 9 cm | 3 / Bag |

BATRIK

REUSABLE PRECISION BRUSH

INSTRUMENT CLEANING

Our Reusable Precision Instrument Cleaning Brush is the ideal brush for the quick removal of difficult debris. It cleans externally as well as deep into instrument jaws. Although the firm and durable bristles feel like stainless steel, there is absolutely no risk of instrument damage or personal injury.

- Bristles have an extra firm density and superior durability
- Polyamide (TYNEX®) bristles
- Steamable thermoplastic ergonomic handle
- Steam / chemical sterilizable
- Latex free



GOOD TO KNOW: ELIMINATES THE NEED FOR THE STAINLESS STEEL TOOTHBRUSH-STYLE BRUSH, WHICH CAN DAMAGE VALUABLE INSTRUMENTS!

PRODUCT SELECTION

| ITEM # | DESCRIPTION | TOTAL LENGTH |
|----------|--|----------------|
| 50-40800 | Reusable Precision External Cleaning Brush | 6.29" / 160 mm |

PACKAGING: Bag of 3

NEW!

TRIFECTA™ BRUSH

INSTRUMENT CLEANING

Introducing the **Trifecta™ Brush**, a versatile **3-in-1 tool** designed to meet all your instrument cleaning needs. This reusable instrument brush is available in two styles: one with **all-nylon** bristles for gentle cleaning and another with nylon bristles on two sides and **stainless steel bristles** on the third, ideal for tackling tougher buildup and stubborn cleaning tasks. This innovative design ensures effective, thorough maintenance for a wide range of instruments.

TRIFECTA™ BRUSH - NYLON

1 NYLON BRISTLES

Gentle, transparent nylon bristles provide an effective yet safe cleaning solution. This side allows you to verify cleanliness easily, without risking damage to delicate surfaces. Ideal for a light, non-invasive clean.

2 PRECISION TIP

Reach those hard-to-clean spots with the precision tip. This small, firm-bristle tip is designed for tight, hard-to-reach areas, such as instrument jaws and box joints. It ensures a deep clean where it's most needed, without the risk of damage or injury to delicate parts.

3 FIRM NYLON BRISTLES

For tougher residues, the firm nylon side delivers the cleaning power you need. With bristles that feel similar to stainless steel but are safe for your instruments, this side easily handles stubborn dirt and buildup.



TRIFECTA™ BRUSH - NYLON & STAINLESS STEEL



PRODUCT SELECTION

| ITEM # | COLOR | DESCRIPTION | LENGTH TOTAL | BRISTLE LENGTH | BRISTLES 1 Ø | BRISTLES 2 Ø | BRISTLES 3 Ø |
|----------|-------|--------------------------------------|----------------|---|-----------------|-----------------|-----------------|
| 28-23012 | Blue | Triple Nylon Brush | 8.66" / 220 mm | 1.20" x 1.20" x 0.47" / 30.5 x 30.5 x 12.0 mm | 0.47" / 12.0 mm | 0.47" / 12.0 mm | 0.315" / 8.0 mm |
| 29-23012 | Blue | Triple Nylon / Stainless Steel Brush | 8.66" / 220 mm | 1.20" x 1.20" x 0.47" / 30.5 x 30.5 x 12.0 mm | 0.47" / 12.0 mm | 0.47" / 12.0 mm | 0.47" / 12.0 mm |

PACKAGING: Bag of 3 - For a bag of 50, add "-50" following the item code.


BATRIK BOTTLE & DRAIN CLEANING BRUSHES

Designed specifically for hospital applications, our Bottle and Drain Cleaning Brushes feature long and flexible bristles to maximize the effectiveness of large-surface cleaning.

- High quality polyamide clear bristles
- Easy to identify debris
- Autoclavable
- Latex free

REUSABLE JAR/BOTTLE BRUSH

STERILIZER DRAIN BRUSH



PRODUCT SELECTION

BLUE DELRIN PLASTIC

| ITEM # | DESCRIPTION | HANDLE LENGTH | BRUSH LENGTH | BRUSH (Ø) | TOTAL LENGTH |
|----------|---------------------------------|---------------|---------------|--------------|-----------------|
| 50-10295 | Reusable Jar / Bottle Brush | 120 mm ± 2 mm | 120 mm ± 5 mm | 7 mm ± 2 mm | 13.19" / 355 mm |
| 50-10224 | Sterilizer Drain Cleaning Brush | N/A | 4" / 10.1 cm | 2" / 5.08 cm | 24" / 609 mm |

PACKAGING: Bag of 3

BATRIK ACETABULAR REAMER BRUSH

Our double-ended cleaning brush is manufactured exclusively for cleaning acetabular bone reamers. The “Pig Tail” rounded style circular head brush is designed to effectively clean inside the reamer, while the rigid straight brush on the reverse side is used for cleaning individual cutting holes. Special fusion of the brush heads to the hand piece along with its high quality antimicrobial medical grade bristles ensures a durable, reusable, and highly effective cleaning tool.



STRAIGHT BRUSH



ROUND BRUSH



PRODUCT SELECTION

| ITEM # | DESCRIPTION | HANDLE LENGTH | STRAIGHT BRUSH | ROUND BRUSH HEAD |
|------------------|----------------------------------|-----------------------------------|-----------------------------------|-------------------------------|
| 50-ACETAB | Acetabular Reamer Cleaning Brush | 4.72" / 120 mm Ø 0.47" / 12 mm | 1.18" / 30 mm Ø 0.3125" / 8 mm | 2" / 50 mm Ø 1.75" / 45 mm |

PACKAGING: Bag of 3

NOTE: Specifications subject to small variances

BATRIK ALL PURPOSE REUSABLE BRUSHES

Plastic Handle Brushes

We offer a wide selection of brushes for instrument cleaning applications. These brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

In-Brush™ Plastic Handle

Available in regular toothbrush style (nylon or stainless) and double-ended brushes.



28-72550



28-69394



29-72550



29-69394



29-70050



50-40911

PRODUCT SELECTION

| ITEM # | DESCRIPTION | LENGTH TOTAL | BRUSH LENGTH | BRUSH HEIGHT |
|-----------------|--|----------------|-----------------------------|------------------------------|
| 28-72550 | Plastic Single End Nylon Reusable Handle Brush | 7.25" / 184 mm | 1.5" / 38 mm | 0.5" / 12.7 mm |
| 29-72550 | Plastic Single End Stainless Steel Reusable Handle Brush | 7.25" / 184 mm | 1.5" / 38 mm | 0.5" / 12.7 mm |
| 29-70050 | Plastic Single End Brass Reusable Handle Brush | 7.25" / 184 mm | 1.5" / 38 mm | 0.5" / 12.7 mm |
| 28-69394 | Plastic Double End Nylon Reusable Handle Brush | 6.9" / 175 mm | 1.38" / 35 mm 1" / 25 mm | 0.39" / 10 mm 0.2" / 5 mm |
| 29-69394 | Plastic Double End Stainless Steel Reusable Handle Brush | 6.9" / 175 mm | 1.38" / 35 mm 1" / 25 mm | 0.39" / 10 mm 0.2" / 5 mm |
| 50-40911 | Hand Scrub Nylon Brush Reusable autoclavable | 4.09" / 104 mm | 3.74" / 95 mm | 1.1" / 28 mm |

PACKAGING: Bag of 3 - For a bag of 50, add "-50" following the item code.

BATRIK

AUTOCLAVE CLEANING BRUSH KITS



ROUND BRUSH

- Features 3 holes for adjustable pole placement
- Adjustable brush placement (parallel or perpendicular) and longer brush length provides easy access to hard-to-reach corners
- Made with FDA approved polyester bristles

HAND BRUSH

- Ergonomic handle for better hand comfort
- Longer, FDA approved polyester bristles
- Easy to clean (no cavity)

PRODUCT SELECTION

REGULAR AUTOCLAVE CLEANING BRUSH KIT

| ITEM # (KIT) | DESCRIPTION | ITEM # (UNIT) | DESCRIPTION | DIMENSIONS |
|--------------|--------------------------------------|---------------|-------------------------------------|--|
| 50-CLAVEREG | Autoclave REGULAR Brush Cleaning Kit | 50-35900 | 1 Hand Brush | Length: 3.5" x 8" / 8.89 cm x 20.32 cm |
| | | 50-55800 | 1 Rounded Corner Brush (Head Brush) | Head size: 5.5 x 10" / 14 x 25.4 cm |
| | | 50-EXTPOLE | 1 Anodized Aluminum Extension Pole | 2-4 ft / 60-122 cm |

SUPERIOR AUTOCLAVE CLEANING BRUSH KIT

| ITEM # (KIT) | DESCRIPTION | ITEM # (UNIT) | DESCRIPTION | DIMENSIONS | CASE DIMENSIONS |
|--------------|---------------------------------------|---------------|-------------------------------------|---------------------------------------|---------------------------------|
| 50-CLAVESUP | Autoclave SUPERIOR Brush Cleaning Kit | 50-35900 | 1 Hand Brush | Length: 3.5" x 8" / 11.5 cm x 20.32cm | |
| | | 50-55800 | 1 Rounded Corner Brush (Head Brush) | Head size: 5.5 x 10" / 14 x 25.4 cm | |
| | | 50-EXTPOLE | 1 Anodized Aluminum Extension Pole | 2-4 ft / 60-122 cm | |
| | | SF-240 | 2 SAVEFACE™ Face Shields | 9" / 23 cm | CASE OF 240 / 26.5" x 21" x 15" |

REPLACEMENT COMPONENTS



HAND BRUSH
50-35900

Packaging: 1



BRUSH HEAD
50-55800

Packaging: 1



EXTENSION POLE
50-EXTPOLE

Packaging: 1



SAVEFACE™
FACE SHIELDS SF-240

Packaging: 100 / Box

BATRIK HYSTEROSONOGRAPHY CATHETER

Hystero-sonography (also known as sonohysterography) can now be performed safely, simply and atraumatically with the gentle infusion of saline into the uterus with the use of our H/S Catheter.


- Available in two sizes: 5FR and 7FR
- Minimally invasive and ideal for office use
- Gentle placement sheath included (decreasing the need for a tenaculum) to minimize patient discomfort
- Flexible and easy-to-use
- Total length of 410mm with a working length of 280 mm
- Latex Free

GOOD TO KNOW: ACHIEVE
IMPROVED IMAGERY
AND DIAGNOSIS WITH
SONOHYSTEROGRAMS

“ We depend upon the reliability, precision and functionality of the Batrik H/S Catheter in our clinical and research operations.

Having performed more than 100,000 sonohysterography procedures, we can unequivocally say that this is the best HSG catheter that we have ever used

”
-Dr. Alex Hartman




PRODUCT SELECTION

| ITEM # | DESCRIPTION | PACKAGING |
|----------------|-------------------|-----------|
| 51-5005 | 5 FR H/S Catheter | 10 / box |
| 51-5007 | 7 FR H/S Catheter | 10 / box |

NOTES





SURGMED GROUP

2165 46th Ave., Lachine, Quebec H8T 2P1
Phone: +1 514-631-7988 | Toll Free: 800-661-5432 | Fax.: +1 514-631-9083
Email: info@surgmed.com | Website: www.batrik.com