




ENDOSCOPY PORTFOLIO

Established in 1997. Proudly Canadian. 
Sold in over 85 countries worldwide.



BATRIK G.O.L.F.F.™ ANTI-FOG SOLUTION

Our G.O.L.F.F.™ Anti-Fog Solution solves the fogging problem often encountered during Laparoscopic surgical procedures.

- Prevents the occurrence of fog on lenses
- Ensures clearer vision and avoids unnecessary interruptions during surgery
- Set includes an easy-squeeze 6 mL bottle and an adhesive radiopaque sponge for secure placement
- OEM (available)
- Sterile
- Latex free



PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING	CASE DIMENSIONS
G.O.L.F.F.	Anti-Fog Solution & Adhesive Sponge set	20 / Box	20" x 16" x 18"


BATRIK ENDOSCOPY BO2™ BITE BLOCKS WITH NASAL OXYGEN DELIVERY

Our BO2™ Bite Blocks significantly enhance patient safety by providing supplemental oxygenation during gastrointestinal endoscopy, preventing hypoxia in sedated patients, and ensuring superior patient comfort.

- Saves time and money by eliminating the need for separate nasal prongs and oxygen masks
- Designed to be free of sharp edges to ensure patient comfort and safety
- Available exclusively in Single-Use
- Optional retention strap is available to ensure a secure fit
- The Bite Blocks and optional retention straps are all latex free, lead free and DEHP free

BO2
(FRONT SIDE)

BO2
(BACK SIDE)



PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING	CASE DIMENSION
BO2	Single-Use	18mm / 54FR, includes strap	100 / Box	10" x 10" x 10"
BO2-N	Single-Use	18mm / 54FR, without strap	100 / Box	10" x 10" x 10"

BATRIK ENDOSCOPY BITESTRONG™ BITE BLOCKS

ADULT & PEDIATRIC BITE BLOCKS

Our Bitestrong™ Bite Blocks provide a safe passage through the patient's mouth during TEE, protecting the patient's teeth and endoscope.

Our molding process ensures a smooth finish with no sharp edges (flashing) that can often irritate a patient's mouth. All Bite Blocks are individually wrapped in sealed, easy to open PE bags and are latex free, lead free, and DEHP free. Additional straps are available, if required

SINGLE USE

- Available without dental retention

REUSABLE

- Designed for repeated usage and will withstand commonly practiced chemical and steam sterilization cycles (134°C / 275°F).
- Available with or without an adjustable retention strap



BS



B1



B3



B4



B6

PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING	CASE DIMENSIONS
B1	Single-Use Pediatric	16 mm / 48FR, includes strap	100 / Box	10" x 10" x 10"
B1-25	Single-Use Pediatric	16 mm / 48FR, includes strap	25 / Box	10" x 7" x 7"
B1-N	Single-Use Pediatric	16 mm / 48FR, without strap	100 / Box	10" x 10" x 10"
B3	Single-Use	20 mm / 60FR, includes strap	100 / Box	10" x 10" x 10"
B3-N	Single-Use	20 mm / 60FR, without strap	100 / Box	10" x 10" x 10"
B3-SA	Single-Use	20 mm / 60FR, with strap <u>attached</u> to one side	100 / Box	10" x 10" x 10"
B3-NDR	Single-Use	20 mm / 60FR, no dental retention rim	100 / Box	10" x 10" x 10"
B3-NDR-N	Single-Use	Same as B3-N, no dental retention rim, without strap.	100 / Box	10" x 10" x 10"
B4	Single-Use	19 mm / 60FR	100 / Box	10" x 10" x 10"
B4T	Single-Use	19 mm / 60FR, includes strap	100 / Box	10" x 10" x 10"
B4T-N	Single-Use	19 mm / 60FR, without strap	100 / Box	10" x 10" x 10"
B6	Reusable	19 mm / 60FR	100 / Box	10" x 10" x 10"
B6T	Reusable	19 mm / 60FR, includes strap	100 / Box	10" x 10" x 10"
B6T-N	Reusable	19 mm / 60FR, without strap	100 / Box	10" x 10" x 10"
B6-25	Reusable	19 mm / 60FR	25 / Box	10" x 7" x 7"
B6T-25	Reusable	19 mm / 60FR, includes strap	25 / Box	10" x 7" x 7"
B6T-N-25	Reusable	19 mm / 60FR, without strap	25 / Box	10" x 7" x 7"
BS	Single-Use	Bite-Block straps 21.7" x 0.79" / 55 cm x 2 cm	100 / Box	10" x 13" x 1.5"

BATRIK COUVRETTE™ VALVE COVERS

Our Couvrette™ Valve Covers eliminate the need for manual cleaning and reprocessing of reusable valves. Our single use Valve Covers effectively reduce cross contamination between patients and also reduce fluid leakage from the parts during procedures.

- Easily and securely fit to all valves
- Reduce the potential for errors
- Single-Use
- Non-Sterile
- Latex free



TIP: REUSABLE VALVE COVERS ARE NOTORIOUSLY DIFFICULT TO CLEAN.

SAVE TIME AND REDUCE RISK WITH COUVRETTE™.

COMING SOON!



PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	PACKAGING
16-0002	Biopsy Valve Cover - Olympus®	Blue	100 / Box
16-0004	Biopsy Valve Cover - Pentax®	Red	100 / Box
16-0022	Biopsy Valve Cover - Olympus® & Fujinon®	Light Blue	100 / Box COMING SOON!
16-0044	Biopsy Valve Cover - Pentax®	Pink	100 / Box COMING SOON!


NOTE: Contact us to verify compatibility with your scope

BATRIK DAYBRUSH™ CHANNEL BRUSHES

REUSABLE

Our Reusable DayBrush™ Channel Brushes can be used for cleaning endoscopes for an entire day. In fact, it's been proven to effectively clean up to 10 scopes before the need to discard.*

- Soft tip ensures easy manipulation within the channel with no risk of damage
- The back-end is sealed and softened to prevent fluid from getting trapped inside the tube. The softening process also prevents damage if the cleaning brush is used in reverse.
- Can be steam sterilized.
- Bristle volume is 20mm in length – nearly double most brushes on the market. A large brush provides more surface contact during the cleaning process.
- Significant per use savings when compared to other reusable or disposable channel cleaning brushes.



() Discard after cleaning 10 scopes or sooner if any brush degradation is noticed. After ten applications, the brush may have the appearance that it can continue to clean effectively, but we discourage prolonged use due to risk of cross contamination. Actual number of uses will vary based on handling and processing.*

PRODUCT SELECTION

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
30-18050-T	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
30-23050-T	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
30-23050-2T	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Double

PACKAGING: 50 / Box

OPTIONAL KIT ADD-ONS

ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only, call for details.


"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.

"FK" to receive both the sponge and valve cleaning brushes.

BATRIK ARC™ CHANNEL BRUSHES SINGLE-USE

BATRIK offers one of the largest selection of Endoscopy Channel Cleaning Brushes in the industry. Our brushes are available in a variety of lengths, materials, and styles to accommodate the widest range of endoscopic devices used in the healthcare industry today.

No matter what your requirements, you'll be able to find the brush that suits your specific needs. If you don't see the brush you require, simply call us. We have over 50 different brushes available and can custom fabricate to your precise specifications.



OPTIONAL KIT ADD-ONS

ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only, call for details.

"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.

"FK" to receive both the sponge and valve cleaning brushes.

SERIES 333 & 311



PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
311-23050	90" / 230 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single
311-23050-2	90" / 230 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Double
311-18050	70" / 180 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single
311-10050	39" / 100 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single

PACKAGING: 50 / BOX

5/6 MM TAPERED BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
333-23056	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.073" / 1.85 mm	Single
333-23056-2	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.073" / 1.85 mm	Double

PACKAGING: 50 / Box

SERIES 300 & 322



PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
300-18050	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Single
300-23050	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Single
300-23050-2	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Double

PACKAGING: 50 / Box

5/6 MM TAPERED BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
322-23056	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.070" / 1.75 mm	Single
322-23056-2	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.070" / 1.75 mm	Double

PACKAGING: 50 / Box

SERIES 344



PRODUCT SELECTION

3 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.55 MM - METALLIC PURPLE DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
344-23030	90" / 230 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single
344-18030	70" / 180 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single
344-10030	39" / 100 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single

PACKAGING: 50 / Box

SERIES 35



PRODUCT SELECTION

2.5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
35-23025	90" / 230 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Single
35-12025	47.2" / 120 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Single

PACKAGING: 50 / Box

2 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
35-23020	90" / 230 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single
35-12020	47.2" / 120 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single
35-18020	70.8" / 180 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single

PACKAGING: 50 / Box

BATRIK VALVE / CONTROL BRUSHES


27-05197



27-57787



40-18511



36-10010



27-15311




PRODUCT SELECTION

ITEM #	DESCRIPTION	LENGTH	BRUSH 1 LENGTH	BRUSH 2 LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)
27-05197	Double Trumpet 5/7 White Color	5.00" / 13 cm	0.20" / 5 mm	1.00" / 25 mm	0.20" / 5 mm	0.28" / 7 mm
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-05197-50						
27-57787	Double. Universal 5/11 Yellow Color	6.00" / 15 cm	1.18" / 30 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-57787-50						
40-18511	Double Universal 5/11 Extra Long Blue Color	8.07" / 20.5 cm	1.57" / 40 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm
PACKAGING: Bag of 3 (reusable and can be steam sterilized). Volume Packaging: Bag of 50 use ITEM # 40-18511-50						
36-10010	Single end large Red Color	4.3" / 11 cm	1.57" / 40 mm	N/A	0.39" / 10 mm	N/A
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 36-10010-50						
27-15311	Single end large Blue Color	4" / 10.3 cm	1.57" / 40 mm	N/A	0.39" / 10 mm	N/A
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-15311-50						


BATRIK V-BRUSHES EXTERNAL CLEANING BRUSHES

Our External Cleaning Brushes for endoscopes / optical instruments provide excellent cleaning capabilities. With their unique U-shape bristles and innovative design, our brushes are the ultimate tool for cleaning the exterior of your endoscopic instruments. Free from sharp or metallic ends that can scratch or cause damage to valuable equipment, both single-use and reusable models are latex free and feature ergonomic handles.



PRODUCT SELECTION

ITEM #	DESCRIPTION	BRUSH LENGTH / (Ø)	HANDLE LENGTH / (Ø)	TOTAL LENGTH	PACKAGING
36-20005	Single-Use External Cleaning Brush	70 mm / 25 mm	80 mm / 3 mm	150 mm	50 / Box
36-20010	Reusable External Cleaning Brush	90 mm / 25 mm	120 mm / 12 mm	210 mm	3 / Bag



SURGMED GROUP

2165 46th Ave., Lachine, Quebec H8T 2P1
Phone: +1 514-631-7988 | Toll Free: 800-661-5432 | Fax.: +1 514-631-9083
Email: info@surgmed.com | Website: www.batrik.com