

# SAVEFACE™ DRAPED SHIELDS

**REDUCE PERSONAL INJURY**  
**SUPERIOR INSTRUMENT DECONTAMINATION PROTECTION**

## SAVEFACE DRAPED - FACE SHIELD

In the Sterile Processing Department (SPD)/Central Sterile Supply Department (CSSD), **personal safety** is paramount.

Our Draped Face Shield provides essential **protection** during high-risk instrument decontamination.

## EFFECTIVE PROTECTION IS NON-NEGOTIABLE

SAVEFACE Draped dynamically reduces harmful exposure to caustic chemicals as well as the risk of cross-contamination. Reduce costly injury to hospital staff.

## PPE CATEGORY 2 CERTIFIED

SAVEFACE DRAPED is tested to protect. Wear it with confidence.



PATENT PENDING



INCLUDES CLOSING SYSTEM TO  
 REDUCE RISK OF FLUID CONTACT  
 FROM BELOW THE SHIELD

## SAVEFACE: KEY FEATURES

### ANSI/AAMI PB70 LEVEL 3 COMPLIANT DRAPE:

- Provides a high level of fluid penetration protection.
- Tuck into clothing or pinch together.

### EXTRA-LONG 9"/23CM WRAP-AROUND SHIELD:

- Ensures full coverage, reducing the risk of exposure from different angles.

### COMFORTABLE:

- The soft, self-shaping foam headpiece molds comfortably to your forehead, allowing for extended wear without discomfort.

### OPTICAL CLARITY WITH ANTI-FOG AND ANTI-GLARE

# DON'T MISS THE OPPORTUNITY TO EXPLORE OUR BROADER RANGE OF PROTECTIVE SHIELDS.

Visit our website [www.surgmed.com](http://www.surgmed.com) or consult our detailed literature to find the shield that best fits your needs.

## VERY IMPORTANT

### IS YOUR PRESENT FACE SHIELD CATEGORY 2 CERTIFIED?

Batrik face shields are independently certified to comply with the PPE Category 2 standard. Many facilities are now requesting compliance to this standard.



SF-240



SF-8240



## PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	TOTAL LENGTH	PACKAGING
SF-240	BATRIK SAVEFACE™ no drape	Shield Length: 9" / 23 cm Shield Width: 13" / 33 cm	9" (23 cm)	240 / Case (8 Boxes of 30 Pcs)
SF-8240	BATRIK SAVEFACE™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23 cm Shield Width: 13" / 33 cm Drape Length: 6" / 15 cm	15" (38 cm)	240 / Case (8 Boxes of 30 Pcs)