

SAVEFACE™ DRAPED SHIELDS

REDUCE PERSONAL INJURY
SUPERIOR INSTRUMENT DECONTAMINATION PROTECTION

SAVEFACE DRAPED-FACE SHIELD

In the Sterile Processing Department (SPD) Central Sterile Supply Department (CSSD), personal safety is paramount.

Our Draped Face Shield provides essential protection during high-risk instrument decontamination.

EFFECTIVE PROTECTION IS NON-NEGOTIABLE

SAVEFACE Draped dynamically reduces harmful exposure to caustic chemicals as well as the risk of cross-contamination. Effective face protection is non-negotiable.

PPE CATEGORY 2 CERTIFIED

SAVEFACE DRAPED is tested to protect. Wear it with confidence.



PATENT PENDING



INCLUDES CLOSING SYSTEM TO REDUCE RISK OF FLUID CONTACT FROM BELOW THE SHIELD

SAVEFACE: KEY FEATURES

ANSI/AAMI PB70 LEVEL 3 COMPLIANT DRAPE:

- Provides a high level of fluid penetration protection.
- Tuck into clothing or pinch together.

EXTRA-LONG 9"/23CM WRAP-AROUND SHIELD:

- Ensures full coverage, reducing the risk of exposure from different angles.

COMFORTABLE:

- The soft, self-shaping foam headpiece molds comfortably to your forehead, allowing for extended wear without discomfort.

OPTICAL CLARITY WITH ANTI-FOG AND ANTI-GLARE

DON'T MISS THE OPPORTUNITY TO EXPLORE OUR BROADER RANGE OF PROTECTIVE SHIELDS.

Visit our website www.surgmed.com or consult our detailed literature to find the shield that best fits your needs.

VERY IMPORTANT

IS YOUR PRESENT FACE SHIELD CATEGORY 2 CERTIFIED?

Batrik face shields are independently certified to comply with the PPE Category 2 standard. Many facilities are now requesting compliance to this standard.

SHIELD LENGTH
9" / 23 CM

DRAPE LENGTH
6" / 15 CM



SF-240



SF-8240



PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	TOTAL LENGTH	PACKAGING	CASE DIMENSIONS
SF-240	BATRIK SAVEFACE™ no drape	Shield Length: 9" / 23cm Shield width: 13' / 33cm	9" (23cm)	240 / Case (8 Boxes of 30 Pcs)	26.5" x 21" x 15"
SF-8240	BATRIK SAVEFACE™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23cm Shield width: 13' / 33cm Shield width: 6" / 15cm	58 cm / 23"	240 / Case (8 Boxes of 30 Pcs)	26.5" x 21" x 15"