


BATRIK VASCULAR SILICONE TIES

AVAILABLE IN 4 SIZES
AND 4 COLORS




PRODUCT SELECTION

ITEM #	DIMENSIONS	COLOR
10-1000	0.075" X 16" (0.2 cm x 40 cm)	Mini-Red
10-1001	0.075" X 16" (0.2 cm x 40 cm)	Mini-Blue
10-1002	0.075" X 16" (0.2 cm x 40 cm)	Mini-Yellow
10-1003	0.075" X 16" (0.2 cm x 40 cm)	Mini-White
10-1010	0.120" x 16" (0.3 cm x 40 cm)	Maxi-Red
10-1011	0.120" x 16" (0.3 cm x 40 cm)	Maxi-Blue
10-1012	0.120" x 16" (0.3 cm x 40 cm)	Maxi-Yellow
10-1013	0.120" x 16" (0.3 cm x 40 cm)	Maxi-White
10-1101	0.021" x 16" (0.5 cm x 40 cm)	Extra-Mini-Blue
10-1111	0.200" x 22" (0.5 cm x 56 cm)	Extra-Maxi-Blue
10-1113	0.200" x 22" (0.5 cm x 56 cm)	Extra-Maxi-White
10-1800	0.075" x 18" / 0.2 cm x 46 cm	Mini-Red
10-1801	0.075" x 18" / 0.2 cm x 46 cm	Mini-Blue
10-1802	0.075" x 18" / 0.2 cm x 46 cm	Mini-Yellow
10-1803	0.075" x 18" / 0.2 cm x 46 cm	Mini-White
10-1810	0.12" x 18" / 0.3 cm x 46 cm	Maxi-Red
10-1811	0.12" x 18" / 0.3 cm x 46 cm	Maxi-Blue
10-1812	0.12" x 18" / 0.3 cm x 46 cm	Maxi-Yellow
10-1813	0.12" x 18" / 0.3 cm x 46 cm	Maxi-White

PACKAGING: 2 ties per sterile pouch – 10 pouches per box

Our Vascular Silicone Ties permit safe and smooth retraction, occlusion and identification of vessels during surgical procedures.

Unlike clamps and retractors which can pinch the vessels flat, or even crack them, the BATRIK elliptical shape allows for the gentle handling of arteries, veins, nerves, tendons, spermatic cords, ureters, or whenever delicate retraction / occlusion is required.



- Manufactured of medical grade silicone
- Do not twist on themselves when unpacked
- Easy to separate, eliminating the need for
- Sterile, radiopaque, and latex free cardboard inserts

PRODUCT SELECTION

ITEM #	DIMENSIONS	COLOR
10-1831	0.04" x 18" / 0.1 cm x 46 cm	Extra-Mini-Blue
10-1841	0.20" x 18" / 0.5 cm x 46 cm	Extra-Maxi-Blue
10-1843	0.20" x 18" / 0.5 cm x 46 cm	Extra-Maxi-White
10-3010	0.075" x 32" (0.2 cm x 81 cm)	Maxi-Red
10-3202	0.075" x 32" (0.2 cm x 81 cm)	Mini-Yellow
10-3200	0.075" x 32" (0.2 cm x 81 cm)	Mini-Red
10-3201	0.075" x 32" (0.2 cm x 81 cm)	Mini-Blue
10-3202	0.075" x 32" (0.2 cm x 81 cm)	Mini-Yellow
10-3203	0.075" x 32" (0.2 cm x 81 cm)	Mini-White
10-3210	0.120" x 32" (0.3 cm x 81 cm)	Maxi-Red
10-3211	0.120" x 32" (0.3 cm x 81 cm)	Maxi-Blue
10-3212	0.120" x 32" (0.3 cm x 81 cm)	Maxi-Yellow
10-3213	0.120" x 32" (0.3 cm x 81 cm)	Maxi-White

PACKAGING: 2 ties per sterile pouch – 10 pouches per box