




PERSONAL PROTECTIVE EQUIPMENT, REUSABLE EQUIPMENT, & CLEANING SUPPLIES



RIGOROUSLY TESTED. FULLY CERTIFIED. LIABILITY INSURED.

Advanced Surgi-Pharm Inc. is proudly Canadian, offering a wide-ranging portfolio of top-rated PPE products, to protect frontline workers (healthcare professionals). From **NIOSH approved N95 respirators, gowns, and face shields** all the way down to **shoe covers and hands-free door openers, your protection as well as that of your patients is a top priority.**

When you purchase your PPE, reusable equipment, and cleaning supplies from Advanced Surgi-Pharm Inc., you can be confident knowing you are purchasing from a medical company with **30 years of experience as an industry leader in high-quality, trusted medical devices.** With over 60 employees dedicated to safety and a **rigorous internal Regulatory/Quality Control team,** all our products are **liability insured** and comply with (or exceed) strict **Health Canada and International standards.**

MASKS / RESPIRATORS

- 4 NIOSH APPROVED N95 MODEL L188 PARTICULATE RESPIRATOR
- 5 NIOSH APPROVED N95 MODEL L288 PARTICULATE RESPIRATOR
- 6 KN95 DISPOSABLE PARTICULATE RESPIRATOR

FACE SHIELDS

- 7 SAVEFACE™ FACE SHIELDS
- 8 BATRIK SHIELD™ FACE SHIELD

PROTECTIVE GOWNS, COVERALLS, AND SHOE COVERS

- 9 LEVEL 2 PROTECTIVE GOWNS (NON-STERILE)
- 10 CHEMOTHERAPY PROTECTIVE LAMINATED STITCHED GOWNS (NON-STERILE)
- 11 LEVEL 3 CHEMOTHERAPY PROTECTIVE LAMINATED ULTRASONICALLY WELDED GOWNS (NON-STERILE)
- 12 LEVEL 3 STERILE SURGICAL LAMINATED ULTRASONICALLY WELDED GOWNS
- 13 LEVEL 4 PROTECTIVE LAMINATED ULTRASONICALLY WELDED & REINFORCED GOWNS (NON-STERILE)
- 14 LEVEL 4 STERILE SURGICAL LAMINATED ULTRASONICALLY WELDED & REINFORCED GOWNS
- 15 DISPOSABLE PROTECTIVE LAMINATED COVERALLS
- 16 SHOOZCOVERS™ SHOE COVERS
- 18 SAVEFACE™ SHOE COVERS

REUSABLE EQUIPMENT

- 19 WTFOG™ REUSABLE ANTI-FOG CLOTHS
- 21 ANTIMICROBIAL IV POLES WITH CuVERRO SHIELD™
- 23 PLEXIGLASS PROTECTION PANELS
- 24 STAINLESS STEEL HANDS-FREE DOOR OPENER

CLEANING SUPPLIES

- 25 BIOESQUE® BOTANICAL DISINFECTANT SOLUTION

NIOSH APPROVED N95 MODEL L188 PARTICULATE RESPIRATOR

Designed for high filtration capacity to protect wearers from airborne particles, this respirator is hypoallergenic and latex free. The flat folded-style adapts well to wearer’s facial features and padded nose piece provides comfort when worn over a prolonged period of time. This N95 respirator has been approved by NIOSH (National Institute for Occupational Safety and Health), the global standard and your only choice for protection you can truly trust.



Performance Standard: The respirator is Certified under 42 CFR 84 and NIOSH approval “TC-84A-6973”

Filtration Efficiency: ≥ 95%

Product Type: Disposable Respirator Mask

Certifications: Filtration efficiency ≥ 95%. Inhalation resistance, exhalation resistance, and Sodium Chloride (NaCl) N95 tested and conform to NIOSH criterion. (NIOSH procedure TEB-APR-STP-0007, RCT-APR-STP-003, and TEB-APR-STP-0059). Health Canada compliant. Approved by Health Canada as an authorized medical device for uses related to COVID-19.

FEATURES

- Conforms highest N95 safety standards, NIOSH N95 (National Institute for Occupational Safety and Health)
- Comfortable fit: flat folded-style adapts well to wearer’s facial features and padded nose piece provides comfort when worn over a prolonged period of time

PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
SD-L188-R0	N95 Respirator Mask - Duckbill Style	20 / Box

NIOSH APPROVED N95 MODEL L288 PARTICULATE RESPIRATOR

Designed for **high filtration capacity** to protect wearers from airborne particles, this N95 respirator is **hypoallergenic, latex free**, and features an adjustable **padded nose piece** for a truly comfortable fit. This N95 respirator **has been approved by NIOSH** (National Institute for Occupational Safety and Health), the global standard and your only choice for protection you can truly trust.



Performance Standard: The respirator is Certified under 42 CFR 84 and NIOSH approval “TC-84A-7228”

Filtration Efficiency: ≥ 95%

Certifications: Filtration efficiency ≥ 95%. Inhalation resistance, exhalation resistance, and Sodium Chloride (NaCl) N95 tested and conform to NIOSH criterion. (NIOSH procedure TEB-APR-STP-0007, RCT-APR-STP-003, and TEB-APR-STP-0059). Health Canada compliant. Approved by Health Canada as an authorized medical device for uses related to COVID-19.

FEATURES

- Conforms to highest N95 safety standards, NIOSH N95 (National Institute for Occupational Safety and Health)
- Comfortable fit: adjustable padded nose piece for comfort when worn over a prolonged period of time

PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
SD-L288-R0	N95 Respirator Mask - Cone Style	20 / Box

BATRIK KN95 DISPOSABLE PARTICULATE RESPIRATOR

The Batrik KN95 Disposable Particulate Respirator is designed to protect wearers from airborne particles. The **latex free** respirator is made with spunbond non-woven and melt blown layers and features an **ear-loop and adjustable nose bridge clip** designed to provide the ultimate in comfort when worn over a prolonged period.



Performance Standard: GB2626 - 2006

Filtration Efficiency: $\geq 95\%$

Certifications: Bacterial filtration efficiency $\geq 95\%$. Inhalation and exhalation resistance tested and conform to NIOSH criterion (NIOSH procedure TEB-APRSTP-0007 & TEB-APR-STP-0003)

FEATURES

- Secured seal: adjustable nose piece for optimize seal. Advanced non-woven fabric makes it easy to breathe.
- Ergonomic Design: durable & stretchable straps with dual point attachments, ensuring reliable protection.

PRODUCT SELECTION


ITEM #	DESCRIPTION	PACKAGING
SF-KN95	KN95 Particulate Respirator	40 / Case

BATRIK SAVEFACE™ FACE SHIELDS

Our Batrik SaveFace™ Face Shields (Visors) are lightweight and allow for good peripheral vision while acting as a physical barrier against splashes. Can be used for extended periods due to comfort-focused design and superior quality materials.

Choose amongst non-draped and draped models:

- Our classic **non-draped** models close the gap between the shield and the chest bone when head is tilted forward for added protection.
- Our **draped** models offer full-length coverage extending from the top of the eyebrow to below the chest, covering the anterior neck area as per CDC recommendations.



HEADPIECE

- Self-shaping foam headpieces mold and cushion forehead for a comfortable and secure fit (per CDC recommendations)
- Shielded air vents maximize airflow for enhanced comfort and prevent fogging

VISOR

- Excellent Distortion-Free Optical Clarity with flat viewing area reduces distortion and glare
- Features anti-fog treatment for clear view
- Wide visors provide peripheral protection to help prevent splashes from going into the eyes
- High Quality PET
 - SF-100/SF-180: 7.5 mil (0.19 mm) thick
 - SF-240/SF-8240: 9.84 mil (0.25 mm) thick

DRAPE

- Our patent-pending drape models come equipped with integrated ANSI/AAMI PB70 Level 3 compliant splash drape barriers for added protection

INCLUDES CLOSING SYSTEM TO REDUCE RISK OF FLUID CONTACT
 PATENT PENDING

ASSEMBLED IN CANADA

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	TOTAL LENGTH	PACKAGING
SF-100	SAVEFACE™ no drape	Shield Length: 9" / 23 cm Shield Width: 13" / 33 cm	9" (23 cm)	100 / Box
SF-180	SAVEFACE™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23 cm Shield Width: 13" / 33 cm Drape Length: 6" / 15 cm	15" (38 cm)	100 / Box
SF-240	BATRIK SAVEFACE™ no drape	Shield Length: 9" / 23 cm Shield Width: 13" / 33 cm	9" (23 cm)	240 / Case (8 Boxes of 30 Pcs)
SF-8240	BATRIK SAVEFACE™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23 cm Shield Width: 13" / 33 cm Drape Length: 6" / 15 cm	15" (38 cm)	240 / Case (8 Boxes of 30 Pcs)

Ver. 2022-01-14

BATRIK SHIELD™ FACE SHIELD

A cost-effective alternative to our premium SAVEFACE™ models, the BATRIK Shield™ provides excellent splash protection. An ideal safety addition to all pandemic kits.

- Self-shaping foam head piece molds comfortably to forehead
- Shields can be used for longer periods due to comfort and improved materials
- Optical shield with anti-fog and anti-glare provide a clearer viewing area
- Soft plastic fasteners adhere the elastic band to the shield for additional comfort and easy use
- Latex free



ASSEMBLED IN CANADA 

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	TOTAL LENGTH	PACKAGING
SF-096	BATRIK Shield™ no drape	Shield Length: 9" / 23 cm*	9" (23 cm)*	96 / Box
SF-896	BATRIK Shield™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23 cm* Drape Length: 7" / 17.8 cm	15.9" / 40.3 cm*	96 / Box

*Shield product tolerances: 8.9" - 9.5" / 22.5 cm - 23.2 cm

BATRIK LEVEL 2 PROTECTIVE GOWNS

LEVEL 2
NON-STERILE
 ULTRASONICALLY WELDED

Our Level 2 Disposable Protective Water Repellant Gowns are designed to protect health professionals and patients from cross contamination.

- Compliant to EN ISO 22610
- Compliant to EN ISO 20811 LIQUID BARRIER PERFORMANCE
- SMS 40 gr/m² soft non-woven spunbond fabric provides protection, breathability, and comfort
- Knitted soft wrist cuffs
- Tie closure
- Ultrasonically welded seams
- Non-Sterile



PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
SF-202-M	M	115 cm / 45"	58 cm / 23"	145 cm / 57"
SF-202-L	L	127 cm / 50"	60 cm / 24"	150 cm / 59"
SF-202-XL	XL	135 cm / 53"	62 cm / 24"	150 cm / 59"
SF-202-XXL	XXL	140 cm / 55"	64 cm / 25"	155 cm / 61"

PACKAGING: 50 GOWNS / BOX

BATRIK

CHEMOTHERAPY PROTECTIVE LAMINATED STITCHED GOWNS

NON-STERILE
STITCHED SEAMS

Our Disposable Chemotherapy Protective Laminated Gowns are impermeable to fluids, offering exceptional protection to healthcare workers handling chemotherapy drugs.

- Fabric Compliant to EN 20811, EN ISO 22610 and ASTM F739-12
- Chemo Resistance to 12 chemotherapy drugs: BiCNU, Adriamycin, Cytosan, Taxol, Toposar, 5-Flu, Thiotepa, Cisplatin, DTIC, Mitomycin C, Oncovin, and Amethopterin
- 56 gr/m² fabric for maximum protection, breathability, and comfort
- Knitted soft wrist cuffs
- Tie closure
- Stitched seams
- Non-Sterile



PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
SF-302-S	S	110 cm / 43"	55 cm / 22"	140 cm / 55"
SF-302-M	M	115 cm / 45"	58 cm / 23"	145 cm / 57"
SF-302-L	L	127 cm / 50"	60 cm / 24"	150 cm / 59"
SF-302-XL	XL	135 cm / 53"	62 cm / 24"	150 cm / 59"
SF-302-XXL	XXL	140 cm / 55"	64 cm / 25"	155 cm / 61"
SF-302-XXXL	XXXL	145 cm / 57"	66 cm / 26"	160 cm / 63"

PACKAGING: 50 GOWNS / BOX

BATRIK

CHEMOTHERAPY PROTECTIVE LAMINATED ULTRASONICALLY WELDED GOWNS

LEVEL 3
NON-STERILE
 ULTRASONICALLY WELDED

Our Level 3 Disposable Chemotherapy Protective Laminated Gowns are impermeable to chemotherapy drugs and other fluids, offering outstanding protection.

- Fabric Compliant to EN 20811 and EN ISO 22610
- Tested and compliant to AAMI PB70 - Level 3 Liquid Barrier Performance
- Tested and compliant to ASTM F739-12
- Chemo Resistance to 12 chemotherapy drugs: BiCNU, Adriamycin, Cytoxan, Taxol, Toposar, 5-Flu, Thiotepa, Cisplatin, DTIC, Mitomycin C, Oncovin, and Amethopterin
- 56 gr/m2 soft non-woven spunbond fabric coated with a polyethylene film provides maximum protection, breathability, and comfort
- Knitted soft wrist cuffs
- Tie closure
- Ultrasonically welded seams
- Non-Sterile



We have been seeing an increasing demand in this smaller yet important market segment

Our products are Product Liability Insured

Our Chemo resistance testing were made at a reputable North American laboratory. Test reports available in English on demand

PRODUCT SELECTION

*ITEM #	**ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
SF-502-S	SF-602-S	S	110 cm / 43"	55 cm / 22"	140 cm / 55"
SF-502-M	SF-602-M	M	115 cm / 45"	58 cm / 23"	145 cm / 57"
SF-502-L	SF-602-L	L	127 cm / 50"	60 cm / 24"	150 cm / 59"
SF-502-XL	SF-602-XL	XL	135 cm / 53"	62 cm / 24"	150 cm / 59"
SF-502-XXL	SF-602-XXL	XXL	140 cm / 55"	64 cm / 25"	155 cm / 61"
SF-502-XXXL	SF-602-XXXL	XXXL	145 cm / 57"	66 cm / 26"	160 cm / 63"

* PACKAGING SF-502: 40 INDIVIDUALLY BAGGED GOWNS / BOX

** PACKAGING SF-602: 50 BULK PACKED GOWNS / BOX

STERILE SURGICAL LAMINATED ULTRASONICALLY WELDED GOWNS

**LEVEL 3
STERILE
ULTRASONICALLY WELDED**

Our Level 3 Disposable Sterile Surgical Laminated Gowns provide high performance with a fluid resistant barrier protection.

- Fabric Compliant to EN 20811 and EN ISO 22610
- Tested and compliant to AAMI PB70 - Level 3 Liquid Barrier Performance
- Tested and compliant to ASTM F739-12
- Chemo Resistance to 12 CHEMOTHERAPY DRUGS: BiCNU, Adriamycin, Cytoxan, Taxol, Toposar, 5-Flu, Thiotepa, Cisplatin, DTIC, Mitomycin C, Oncovin, and Amethopterin
- 56 gr/m² fabric for maximum protection, breathability, and comfort
- Knitted soft wrist cuffs
- Tie closure
- Ultrasonically welded seams
- Sterile and aseptically packed

Individually packed, each gown comes in a set designed to be used during surgery to reduce contamination and infection between the patient and medical team.



Each set includes:

- 1 x Sterile Surgical Laminated Gown - Level 3 (Size S - XXXL)
- 1 x Blue Sterilization Wrap (60 cm x 60 cm / 16" x 16")
- 1 x White Cellulose Hand Towel (30 cm x 40 cm / 12" x 16")

PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
41.001.01	S	110 cm / 43"	55 cm / 22"	140 cm / 55"
41.002.01	M	115 cm / 45"	58 cm / 23"	145 cm / 57"
41.003.01	L	127 cm / 50"	60 cm / 24"	150 cm / 59"
41.004.01	XL	135 cm / 53"	62 cm / 24"	150 cm / 59"
41.005.01	XXL	140 cm / 55"	64 cm / 25"	155 cm / 61"
41.006.01	XXXL	145 cm / 57"	66 cm / 26"	160 cm / 63"

PACKAGING: 40 INDIVIDUALLY BAGGED SETS / BOX

BATRIK

LEVEL 4 PROTECTIVE LAMINATED ULTRASONICALLY WELDED & REINFORCED GOWNS

LEVEL 4
NON-STERILE
 ULTRASONICALLY WELDED
 & REINFORCED SEAMS

Our **Level 4 Disposable Protective Laminated Gowns** provide high performance for use in high risk situations and/or when infectious diseases are suspected.

- Fabric Compliant to **EN 20811** and **EN ISO 22610**
- Tested and compliant to **AAMI PB70 - Level 4 Liquid Barrier Performance and blood and viral penetration resistance**
- 56 gr/m² soft non-woven spunbond fabric coated with a polyethylene film provides maximum **protection, breathability, and comfort**
- Tested and compliant to **ASTM F739-12**
- Knitted soft wrist cuffs
- Tie closure
- **Ultrasonically welded & reinforced seams**
- **Non-Sterile**



**INSIDE VIEW OF
ULTRASONICALLY WELDED &
REINFORCED SEAMS**

PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
SF-402-S	S	110 cm / 43"	55 cm / 22"	140 cm / 55"
SF-402-M	M	115 cm / 45"	58 cm / 23"	145 cm / 57"
SF-402-L	L	127 cm / 50"	60 cm / 24"	150 cm / 59"
SF-402-XL	XL	135 cm / 53"	62 cm / 24"	150 cm / 59"
SF-402-XXL	XXL	140 cm / 55"	64 cm / 25"	155 cm / 61"
SF-402-XXXL	XXXL	145 cm / 57"	66 cm / 26"	160 cm / 63"

PACKAGING: 50 GOWNS / BOX

STERILE SURGICAL LAMINATED ULTRASONICALLY WELDED & REINFORCED GOWNS

**LEVEL 4
STERILE**
ULTRASONICALLY WELDED
& REINFORCED SEAMS

Our Level 4 Disposable Sterile Surgical Laminated Gowns provide high performance with a fluid resistant barrier protection.

- Fabric Compliant to EN 20811 and EN ISO 22610
- Tested and compliant to AAMI PB70 - Level 4 Liquid Barrier Performance and blood and viral penetration resistance
- 56 gr/m² soft non-woven spunbond fabric coated with a polyethylene film provides maximum protection, breathability, and comfort
- Tested and compliant to ASTM F739-12
- Knitted soft wrist cuffs
- Tie closure
- Ultrasonically welded & reinforced seams
- Sterile and aseptically packed

Individually packed, each gown comes in a set designed to be used during surgery to reduce contamination and infection between the patient and medical team.

Each set includes:

- 1 x Sterile Surgical Laminated Gown - Level 4 (Size S - XXXL)
- 1 x Blue Sterilization Wrap (60 cm x 60 cm / 16" x 16")
- 1 x White Cellulose Hand Towel (30 cm x 40 cm / 12" x 16")



INSIDE VIEW OF
ULTRASONICALLY
WELDED &
REINFORCED SEAMS

PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
46.001.01	S	110 cm / 43"	55 cm / 22"	140 cm / 55"
46.002.01	M	115 cm / 45"	58 cm / 23"	145 cm / 57"
46.003.01	L	127 cm / 50"	60 cm / 24"	150 cm / 59"
46.004.01	XL	135 cm / 53"	62 cm / 24"	150 cm / 59"
46.005.01	XXL	140 cm / 55"	64 cm / 25"	155 cm / 61"
46.006.01	XXXL	145 cm / 57"	66 cm / 26"	160 cm / 63"

PACKAGING: 40 INDIVIDUALLY BAGGED SETS / BOX

DISPOSABLE PROTECTIVE LAMINATED COVERALLS

Our Disposable Protective Laminated Coveralls provide comfort and outstanding protection for high risk situations (Type 5-B & 6-B)

NON-STERILE
OVERLOCK STITCHED
SEAMS TYPE 5-B & 6-B

- Compliant to EN 14126, EN ISO 13982-1 (Type 5) and EN 13034 (Type 6)
- Compliant to ISO 16603 and ISO 16604 for blood and viral penetration
- Tested for puncture resistance (EN 863) and tested against chemicals (EN 14325)
- Biflex 55 gr/m² soft non-woven spunbound material coated with a polyethylene film provides maximum protection, breathability, and comfort
- Inside overlock stitched seams
- Zipper with storm flap
- Elastic hood, cuffs, ankles, and waist
- Non-Sterile



PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	CHEST WIDTH
95.00S.02	S	162 cm - 170 cm / 64" - 67"	84 cm - 92 cm / 33" - 36"
95.00M.02	M	168 cm - 176 cm / 66" - 69"	92 cm - 100 cm / 36" - 39"
95.00L.02	L	174 cm - 182 cm / 68.5" - 72"	100 cm - 108 cm / 39" - 42.5"
95.00XL.02	XL	180 cm - 188 cm / 73" - 74"	108 cm - 116 cm / 42.5" - 46"
95.00XXL.02	XXL	186 cm - 194 cm / 73" - 76"	116 cm - 124 cm / 46" - 49"

PACKAGING: 40 INDIVIDUALLY POLYBAGGED COVERALLS / BOX

BATRIK

SHOOZCOVERS™

Batrik ShoozCovers™ are manufactured in Canada from high quality 35 g/m² polypropylene fabric and 100% latex free elastics for performance and durability you can rely on.


Our anti-skid tread reduces the risk of slipping and our smart design delivers easy on-and-off convenience and great fit.

WHY YOU'LL LOVE THEM:

- **Low linting non-woven 35 g/m² polypropylene fabric** protects footwear and keeps floors clean
- Particle count tested and complies with IEST-RP-CC003.4
- Shooz-Repellent™ technology delivers **basic incidental fluid contact repellency**¹
- **Triple ultrasonic welds** increase tear strength and reduce linting
- High quality **made in the USA latex-free elastics** are triple reinforced to deliver a strong and secure fit
- **Generous sizing** and smart design for better fit and easy on-and-off convenience
- **Anti-skid tread on soles**
- **Truly 100% latex and silicone free shoe covers**²
- **Non-Sterile, Single Use**
- **Long-lasting to reduce waste**



MADE IN CANADA 



PERFECT FOR:

- Protecting a variety of footwear from shoes to clogs to boots
- Use in medical establishments, pharmaceutical labs, electronic manufacturing, GMP sites, and anywhere optimum hygiene is required
- Transitioning from “unclean” to “clean” areas

¹Shooz-Repellent™ technology provides basic incidental fluid contact repellency. It is not waterproof nor water resistant.

²Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product.

TWO VERSATILE SIZES:

UNIVERSAL


UN

Up to shoe size 11 (44.5 EU)
~ 42 cm x 16.5 cm / 16.5" x 6.5"³

EXTRA LARGE

XL

Shoe sizes 11-15 (44.5-49.5 EU)
~ 47 cm x 19 cm / 18.5" x 7.5"³



SHOOZCovers™ are made longer (42 & 47cm) for accurate shoe size selection.

DONNING IS AS EASY AS 1-2-3:



MADE IN CANADA

- Support Canadian manufacturing
- Eliminate supply chain issues
- Reduce shipping costs

PRODUCT SELECTION

ITEM #	DESCRIPTION	SIZE	PACKAGING
SC-1216-PDB	ShoozCovers™ Disposable Shoe Covers Blue 35 g/m ² fabric	Universal Up to shoe size 11 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	100 / Dispenser Box 6 Dispenser Boxes / Case
SC-1226-PDB	ShoozCovers™ Disposable Shoe Covers Blue 35 g/m ² fabric	XL Shoe size 11 - 15 ~ 47 cm x 19 cm / 18.5" x 7.5"	100 / Dispenser Box 6 Dispenser Boxes / Case

BATRIK

SAVEFACE™ SHOE COVERS

Our SaveFace™ superior medical shoe covers are durable, safe, and easy to use. Made of polypropylene fabric with robust welds, our shoe covers offer the ultimate in durability. The triple-elastic run and anti-skid print allows for a smooth fit around the foot and reduces the risk of slipping.

- Particle count tested and complies with IEST-RP-CC003.4
- 100% latex and silicone free



MADE IN CANADA 

PRODUCT SELECTION

ITEM #	SIZE	PACKAGING
SF-SC02UV	Universal Up to shoe size 11	Regular Packaging 200 / dispenser box - 4 dispenser boxes / case
SF-SC02UV-8	Universal Up to shoe size 11	Volume Packaging 100 / bag - 8 bags / case
SF-SC02XL	Extra-Large Size 11 - 15	Regular Packaging 200 / dispenser box - 4 dispenser boxes / case
SF-SC02XL-8	Extra-Large Size 11 - 15	Volume Packaging 100 / bag - 8 bags / case

BATRIK
WTFOG™ by G.O.L.F.F.™
REUSABLE ANTI-FOG CLOTHS

“WHAT THE FOG?!”

We’ve all experienced it.

Eyewear fogging up while wearing a mask is not only unpleasant and highly frustrating, but it can be dangerous, too.

INTRODUCING THE FOG FIX:

Our WTFog™ reusable anti-fog cloths prevent fogging on all types of eyewear. Ideal for use when wearing a mask or face shield, WTFog™ takes you from “What The Fog?!” to confidently ready for your next patient in a snap.

WHY YOU’LL LOVE IT:

- Eliminates fogging on all glass, polycarbonate, and plastic eyewear lenses
- Each reusable nano-infused cloth lasts up to 600 wipes or 6 months
- Microfiber cloths are infused with our proprietary heavy duty anti-fog treatment
- Quick and easy application - just wipe and wear
- Lint-free, scratch-free, no greasy smearing
- Long lasting protective coating provides effective protection against fog, humidity, sweat, and fast temperature changes.
- 100% Natural Titanium Dioxide (TiO2) Nano Formula developed in Germany meets strict REACH standards
- Non-Sterile



“To be able to complete a whole appointment without adjusting my eye protection makes maintaining COVID PPE protocols that much easier.”

DR. MOLLY FRIEDMAN

HOLLAND BLOORVIEW KIDS REHABILITATION HOSPITAL
 TORONTO, ON

PERFECT WHEN:

- You wear glasses/goggles and a mask/face shield
- Your profession regularly exposes you to fogging and you need crystal clear vision without the risk or hassle of fogging up



DIRECTIONS FOR USE:

Open resealable pouch and remove the WTFOG™ reusable anti-fog cloth. Wipe cloth on clean surface (glasses/goggle lens) being careful not to remove the protective coating delivered by the cloth. There is no wait time required before wearing your protected item. Enjoy your fog-free experience!

CLOTH STORAGE & CARE:

After use, keep cloth in the resealable pouch until next use, ensuring the seal is fully closed. If cloth becomes too dry, add up to 1 ml of water to pouch to reactivate cloth. For best results, store at ambient temperature/humidity. Dispose cloth after 600 uses or 6 months.

External surface use only.
Keep out of reach of children.
Do not wash

PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
WTFOG-30	30 WTFOG™ in a dispenser box measuring: 3 7/16" x 5 2/16" x 6 2/16" (8.7 cm x 13 cm x 15.4 cm)	30 WTFOG™ / dispenser box
WTFOG-100	100 WTFOG™ - bulk packed	100 WTFOG™ / box
WTFOG-1200	1200 WTFOG™ - bulk packed	1200 WTFOG™ / case (12 boxes of 100 bulk packed)

1 WTFOG™ Cloth per Resealable Pouch
WTFOG™ Cloth: 5.5" x 5.5" (13.9 cm x 13.9 cm)
WTFOG™ Resealable Pouch: 3.5" x 5" (8.9 cm x 12.7 cm)

IMPERIAL SURGICAL LTD.

ANTIMICROBIAL IV POLES WITH CUVERRO SHIELD™


Frequently touched surfaces like IV poles have long been identified as reservoirs for the spread of pathogenic microbes.¹ These can easily contaminate the hands and equipment of healthcare professionals who in turn can transmit these pathogens to patients during routine care leading to Healthcare-associated infections.

Our Antimicrobial IV Poles protected with CuVerro Shield™ by Aereus Technologies* offers you a solution to this problem.

CuVerro Shield™ is a unique coating of naturally antimicrobial copper that is applied to our stainless steel poles, ram horns, and u-hooks to deliver outstanding protection to your highest touch zones.

The result? CuVerro® kills more than 99.9% of bacteria** within two hours and continues to kill 99% of bacteria** even after repeated contamination.

A US government funded study conducted by researchers at the US National Institutes of Health and the Centers for Disease Control and Protection found that the HCoV-19 virus remained viable for up to 2 to 3 days on plastic and stainless steel surfaces versus up to 4 hours on copper.²



Viability of SARS-CoV-1 and SARS-CoV-2 in Aerosols and on Various Surfaces.²

* CuVerro Shield™ by Aereus Technologies is a Health Canada Pest Management Regulatory Agency (PMRA) registered process.

**Laboratory testing shows that, when cleaned regularly, CuVerro Surfaces kill greater than 99.9% of the following bacteria within 2 hours of exposure: Methicillin-Resistant Staphylococcus aureus, Staphylococcus aureus, Enterobacter aerogenes, Pseudomonas aeruginosa, E.coli O157:H7, and Vancomycin-Resistant Enterococcusfaecalis (VRE). The use of CuVerro® bactericidal copper products is a supplement to and not a substitute for standard infection control practices; users must continue to follow all current infection control practices, including those practices related to cleaning and disinfection of environmental surfaces. This surface has been shown to reduce microbial contamination, but it does not necessarily prevent cross contamination.

¹Oie S, Hosokawa I, Kamiya A. 2002. Contamination of room door handles by methicillin sensitive/methicillin-resistant Staphylococcus aureus. J. Hosp. Infect. 2002;51:140-143.

²Doremalen NV, Bushmaker T, Morris DH, et al. Aerosol and surface stability of HCoV-19 (SARS-CoV-2) compared to SARS-CoV-1. N Engl J Med 2020.

BASES & CASTERS

OR-36770-CUV	CuVERRO SHIELD™ IV POLE, CENTRE COLLAR AND KNOB with 23" cream aluminium base (V-400016) and 2" casters
OR-36760-25-CUV	CuVERRO SHIELD™ IV POLE, CENTRE COLLAR AND KNOB with 25" stainless steel base (V-400013-25) and 2" casters
OR-36760-253-CUV	with 25" stainless steel base (V-400013-25) and 3" casters
OR-36771-25-CUV	CuVERRO SHIELD™ IV POLE, CENTRE COLLAR AND KNOB with 25" stainless steel base (V-400013-25), 2" casters and weight kit, of 12.5 lbs. (V-50004)
OR-36771-253-CUV	with 25" stainless steel base (V-400013-25), 3" casters and weight kit, of 12.5 lbs. (V-50004)

SMOOTH-ROLLING
NON-MARKING
DURABLE CASTERS




OPTIONAL WEIGHT KIT
AVAILABLE FOR INCREASED
STABILITY

HOOKS & RAM HORNS

OR-36499-CUV	CuVERRO SHIELD™ STAINLESS STEEL RAM HORNS AND ENDCAP Set of 2
OR-36500-CUV	CuVERRO SHIELD™ STAINLESS STEEL RAM HORNS AND ENDCAP Set of 4
OR-36501-CUV	CuVERRO SHIELD™ STAINLESS STEEL U-HOOKS AND ENDCAP Set of 2
OR-36502-CUV	CuVERRO SHIELD™ STAINLESS STEEL U-HOOKS AND ENDCAP Set of 4
OR-36504-CUV	CuVERRO SHIELD™ STAINLESS STEEL RAM HORNS (2), STAINLESS STEEL U-HOOK (2) AND ENDCAP Set of 4
V-50002-CUV	CuVERRO SHIELD™ STAINLESS STEEL URINARY DRAIN HOOK AND COLLAR WITH KNOB Set of 1

HEAVY-DUTY RAM
HORNS, U-HOOKS OR A
COMBINATION OF BOTH
AVAILABLE



POLES & HANDLES

V-CUV-50003-1	CuVERRO SHIELD™ ROUND STEEL HANDLEW
V-CUV-50006-1	CuVERRO SHIELD™ STAINLESS STEEL TRIANGULAR HANDLE
V-CUV-50007	CuVERRO SHIELD™ SEMI-ROUND STEEL HANDLE

LARGE PRONG KNOB FOR
HEIGHT ADJUSTMENT, EASY
TO TURN EVEN WITH WET
HANDS



GOOD TO KNOW: In the unlikely event your CuVerro coated IV Pole should become deeply scratched, a study conducted by the University of Waterloo demonstrated that antimicrobial activity of the Aereus coating is maintained even when completely removed up to a width of between 3 and 5 mm, under the test conditions used in this study.³

³ Tanvir S, Anderson WA. 2020. Antimicrobial Activity of "Deeply Scratched" Aereus-Coated Surfaces. Department of Medical Engineering. University of Waterloo. 2020.


IMPERIAL SURGICAL LTD.

PLEXIGLASS PROTECTION PANELS

Our Plexiglass Protection Panels are made of 100% acrylic clear plastic to safeguard face-to-face interaction. They are available in a variety of ready-made standard sizes plus we can create custom sizes/designs to accommodate your specific needs. Our simple and ergonomic design provides the necessary protection for your patients, customers, and community without an imposing presence.

FEATURES

- Multiple dimensions and opening sizes available
- 100% acrylic clear plastic
- Easy to install and move as needed
- Easy to disinfect



PRODUCT SELECTION

ITEM #	DIMENSIONS	ACRYLIC THICKNESS	OPENING SIZE
PG-2331	23.5" W x 31.5" H	3/16"	12" W x 4" H
PG-2336	23.5" W x 36" H	1/4"	8" W x 2" H
PG-2432	24" W x 32" H	3/16"	12" W x 4" H
PG-2436	24" W x 36" H	3/16"	12" W x 4" H
PG-3531	35.5" W x 31.5" H	3/16"	12" W x 4" H
PG-3632	36" W x 32" H	3/16"	5" W x 10" H
PG-3636	36" W x 36" H	3/16"	5" W x 10" H
PG-3735	37.5" W x 35.5" H	1/4"	8" W x 2" H
PG-4831	48" W x 31.5" H	3/16"	12" W x 4" H


IMPERIAL SURGICAL LTD.

STAINLESS STEEL HANDS-FREE DOOR OPENER

The Imperial Surgical Ltd. Hands-Free Door Opener is a 100% stainless steel, foot operated door opener that works on any latch-less door.

Easy to install and easy to clean, the no-touch opener helps minimize the spread of germs by reducing the need to touch door handles where bacteria and viruses may live.

The hands-free door opener comes with an eye catching sticker to put on the door to encourage people to use it. We also provide the hardware so you have everything you need for installation.



PRODUCT SELECTION

ITEM #	DIMENSIONS	PACKAGING
11656-R0	5" Length x 1.5" Height x 0.14" Thickness (127 mm Length x 38 mm Height x 3.5 mm Thickness)	1 Unit

BIOESQUE® BOTANICAL DISINFECTANT SOLUTION

Capable of cleaning, sanitizing, disinfecting, and deodorizing in a single step, Bioesque's ready to use Botanical Disinfectant Solution is the natural solution for a wide range of cleaning challenges. Free of bleach and other harsh chemicals, it is safe for everyday use. Powered by the patented Thymox technology, it utilizes Thymol, a cutting-edge antimicrobial agent present as a component of Thyme Oil, to efficiently and effectively disinfect.

FEATURES

- Kills 99.9% of Bacteria, Viruses*, Fungi & Molds and eliminates 99.9% of most Allergens*** upon contact
- Bactericidal, Virucidal*, Tuberculocidal, and Fungicidal**
- Kills Norovirus in just 4 Minutes
- No Rinse Required, even on food contact surfaces
- Sanitizes Soft Surfaces in just 2 minutes
- Non-Abrasive, Non-Corrosive, and Non-Flammable
- Free from chlorine, phosphates, phthalates, propylparaben, butylparaben, formaldehyde, formaldehyde-donors, and NPE's
- Pleasant lemongrass grapefruit scent effectively controls strong, unpleasant odors



PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
BIO-UN946-R0	Bioesque botanical disinfectant solution	1 x 946 ml Bottle
BIO-CS946-R0	Bioesque botanical disinfectant solution	12 x 946 ml Bottles
BIO-UNGAL-R0	Bioesque botanical disinfectant solution	1 x Gallon (4/ Case)

*Virucidal according to the ASTM Standard Test Method for Efficacy of Virucidal Agents

** Fungicidal according to the AOAC Germicidal Spray Method

*** Removes the following allergens: dust mite matter, cat dander, dog dander, mouse dander, rat dander, guinea pig dander, cockroach matter, peanut, birch pollen, timothy grass pollen, rag weed pollen.



2165 46th Ave., Lachine, QC, Canada, H8T 2P1
Phone: +1.514.631.7988 | Toll Free: 800.661.5432 | Fax: +1.514.631.9083
Email: info@surgmed.com | Website: www.advancedsurgipharm.com