



Batrik's Essentials

for superior sterilization!



PRIORITY
Room # _____
Time Needed _____
Initials _____

BATRIK

TAG-IT™

Batrik Tag-It™ is our new line of information tags developed to help medical professionals identify and document a wide range of instructions. Whether you need to identify instruments, trays or carts, both before or after sterilization, Tag-It™ has got you covered. These durable and tear-resistant tags are perfect for indicating special handling instructions.

WHY YOU’LL LOVE THEM:

- The paper, tie wraps, and smudge-proof ink are designed to withstand high temperatures of up to 273.2°F (134°C) and are compatible with Gas plasma sterilization (STERRAD®100 NX Sterilization System).
- The synthetic paper is durable and will not wrinkle or absorb water. Made from recyclable plastic promoting environmental sustainability.
- Robust and easy to use tie wraps designed to securely attach to a wide range of trays, instruments, and other medical devices.
- Available in over 8 languages.

KEY PRODUCT FEATURES:

- **High Quality Permanent Marker**
- **Steam or Plasma indicator Dots** provide confirmation of processing
- **Eco-friendly, Waterproof & Tearproof Tags**
- **Large Quantities:** Each bundle of Tag-It™ comes with 100 versatile tags, to quickly identify instruments that require attention.
- **Available in multiple linguistic choices**



TI-3503

All tags available with indicator dots
Provides visual confirmation that the tag was processed



PRODUCT SELECTION:

ITEM #	COLOR	DESCRIPTION	SIZE	PACKAGING
TI-3502	Blue	VENDOR	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3503	Green	LOANER	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3504	Red	PRIORITY	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3505	Yellow	RE-RUN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3506	Pink	MISSING	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3507	Orange	BROKEN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-3509	Black	BLANK	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35B1	White	BLANK	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35B4	Red	BROKEN	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35C22	Light Blue	VENDOR CHECK IN + OUT	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35L3	Green	LOANER CHECK IN + OUT	3" x 5" / 8 cm x 13cm	100 / PKG
TI-35T2	Blue	QUICK TURNOVER	3" x 5" / 8 cm x 13cm	100 / PKG
TI-PM09	Black	FINE PERMANENT MARKERS	1.5" Ø x 5.5" / 3.81cm Ø x 12.7 cm	36 / Bag

All information tags are single-use. **PACKAGING: 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.**

PROCESSING INDICATOR DOTS: Add suffix "ST" to end of code for steam dots.

PROCESSING INDICATOR DOTS: Add suffix "PL" to end of code for plasma dots.

PRODUCT SELECTION:



<p>VENDOR / BLUE</p> <p>ITEM # SIZE</p> <hr/> <p>TI-3502 3" x 5" (8 cm x 13cm)</p>			<p>BLANK / WHITE</p> <p>ITEM # SIZE</p> <hr/> <p>TI-35B1 3" x 5" (8 cm x 13cm)</p>
<p>LOANER / GREEN</p> <p>ITEM # SIZE</p> <hr/> <p>TI-3503 3" x 5" (8 cm x 13cm)</p>			<p>BROKEN / RED</p> <p>ITEM # SIZE</p> <hr/> <p>TI-35B4 3" x 5" (8 cm x 13cm)</p>
<p>PRIORITY / RED</p> <p>ITEM # SIZE</p> <hr/> <p>TI-3504 3" x 5" (8 cm x 13cm)</p>			<p>VENDOR IN + OUT / LIGHT BLUE</p> <p>ITEM #</p> <hr/> <p>TI-35C22</p> <p>SIZE</p> <hr/> <p>5" x 3" (13 cm x 8 cm)</p>
<p>RE-RUN / YELLOW</p> <p>ITEM # SIZE</p> <hr/> <p>TI-3505 3" x 5" (8 cm x 13cm)</p>			<p>LOANER IN + OUT / GREEN</p> <p>ITEM #</p> <hr/> <p>TI-35L3</p> <p>SIZE</p> <hr/> <p>5" x 3" (13 cm x 8 cm)</p>
<p>MISSING / PINK</p> <p>ITEM # SIZE</p> <hr/> <p>TI-3506 3" x 5" (8 cm x 13cm)</p>			<p>BROKEN / ORANGE</p> <p>ITEM # SIZE</p> <hr/> <p>TI-3507 3" x 5" (8 cm x 13cm)</p>
<p>BLANK / BLACK</p> <p>ITEM # SIZE</p> <hr/> <p>TI-3509 3" x 5" (8 cm x 13cm)</p>			<p>QUICK TURNOVER/ BLUE</p> <p>ITEM # SIZE</p> <hr/> <p>TI-35T2 3" x 5" (8 cm x 13cm)</p>

All information tags are single-use. PACKAGING: 1 pack of 100 tags, 100 tie wraps, and 1 permanent marker.


BATRIK SILICONE PIN MAT

Our Silicone Pin Mat cushions, protects, and organizes delicate surgical instruments for sterilization. The elevated nubs provide greater air-flow ventilation and drying, significantly reducing residual moisture collection.

- Reusable
- Autoclavable
- Latex free



FRONT SIDE
Pins keep instruments securely in place



BACK SIDE
Numerous holes and elevated nubs provide superior venting and reduce water collection

PRODUCT SELECTION

ITEM #	DIMENSIONS	PACKAGING
13-2001	13" x 18" (33 cm x 45.72 cm) 11.97 mm high from base of mat to end of pin 15.53 mm high from base of elevated nub to end of pin	1 unit
13-2005	13" x 18" (33 cm x 45.72 cm) 11.97 mm high from base of mat to end of pin 15.53 mm high from base of elevated nub to end of pin	5 units

GOOD TO KNOW:
CUSTOM SIZES
ARE AVAILABLE
UPON REQUEST

BATRIK CORNER POCKETS™


With three models to choose from, our Corner Pockets™ portfolio is a truly comprehensive line of products designed to prevent tearing associated with sharp corners or feet on surgical trays.



18-CP02



18-CP07




18-CP06

COLOUR	BLUE	ORANGE	PINK
USAGE	SINGLE-USE	REUSABLE <i>withstands repeated sterilization up to 134°C</i>	REUSABLE <i>withstands repeated sterilization up to 134°C</i>
DIMENSIONS	2.8" x 4" (7.3 cm x 10.3 cm)	2.8" x 4" (7.3 cm x 10.3 cm)	3.4" x 4.7" (8.5 cm x 12 cm)
LATEX-FREE	✓	✓	✓
NO MOISTURE RETENTION	✓	✓	✓
FABRICATION	Flexible, silicone-like plastic	Flexible, silicone-like plastic	Our most flexible plastic for easiest placement and most secure protection

SOFT POCKETS™

Our Soft Pockets™ Foam Corner Protectors are made from high density foam to effectively cushion and protect the corners of surgical trays to prevent tearing. They are especially useful with heavier instrument sets.



PRODUCT SELECTION




ITEM #	DESCRIPTION	PACKAGING
18-CP02**	Corner Pockets™ Single-Use, Blue	240 / Bag
18-EVAL2**	Corner Pockets™ Single-Use, Blue	24 / Bag
18-CP07	Corner Pockets™ Reusable, Orange	80 / Bag
18-CORNER	Corner Pockets™ Reusable, Orange	320 / Box (4 Bags of 80 pcs)
18-EVAL7	Corner Pockets™ Reusable, Orange	24 / Bag
18-CP06	Corner Pockets™ Reusable, Pink	80 / Bag
18-CPF-02	Soft Pockets™ Single-Use Foam Corner Protectors	80 / Box (24 Boxes / Pallet)

** Although the single use Corner Pockets™ appear similar to the reusable models they are blended differently and have only been validated as single-use. Re-using them may compromise sterility.

BATRIK INSTRUMENT PROTECTORS

Our Instrument Protectors protect the tips and edges of instruments to prevent damage during handling and sterilization. Available in a wide array of color coded sizes for almost every application, including osteotome, retractor, double, and perforated.

- Steam and gas permeable
- Perforated protectors reduce vacuum seal and moisture build-up
- Validated for pre-vacuum sterilization*
- Multiple color-coded sizes
- Latex Free
- Non-Sterile



*Independent laboratory tested to validate pre-vacuum steam sterilization efficacy in a healthcare setting

PRODUCT SELECTION

REGULAR & PERFORATED / CLEAR



RC-07



PC-07

ITEM # REGULAR	ITEM # PERFORATED	COLOR	DIMENSIONS
RC-01	PC-01	Clear	0.062" x 0.75" (1.57 mm x 19 mm)
RC-02	PC-02	Clear	0.075" x 0.75" (2 mm x 19 mm)
RC-03	PC-03	Clear	0.109" x 0.75" (2.77 mm x 19 mm)
RC-04	PC-04	Clear	0.125" x 1" (3.18 mm x 25.4 mm)
RC-05	PC-05	Clear	0.187" x 1" (4.75 mm x 25.4 mm)
RC-06	PC-06	Clear	0.406" x 0.75" (10.3 mm x 19 mm)
RC-07	PC-07	Clear	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
RC-08	PC-08	Clear	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
RC-09	PC-09	Clear	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
RC-10	PC-10	Clear	Assorted

PACKAGING: Bag of 100

REGULAR & PERFORATED / TINTED



RT-07



PT-07

RT-01	PT-01	Tinted White	0.062" x 0.75" (1.57 mm x 19 mm)
RT-02	PT-02	Tinted Blue	0.075" x 0.75" (2 mm x 19 mm)
RT-03	PT-03	Tinted Green	0.109" x 0.75" (2.77 mm x 19 mm)
RT-04	PT-04	Tinted Red	0.125" x 1" (3.18 mm x 25.4 mm)
RT-05	PT-05	Tinted Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
RT-06	PT-06	Tinted Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
RT-07	PT-07	Tinted Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
RT-08	PT-08	Tinted Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
RT-09	PT-09	Tinted Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
RT-10	PT-10	Tinted Assorted	Assorted

PACKAGING: Bag of 100

REGULAR & PERFORATED / SOLID



RS-07




PS-07

RS-01	PS-01	Solid White	0.062" x 0.75" (1.57 mm x 19 mm)
RS-02	PS-02	Solid Blue	0.075" x 0.75" (2 mm x 19 mm)
RS-03	PS-03	Solid Green	0.109" x 0.75" (2.77 mm x 19 mm)
RS-04	PS-04	Solid Red	0.125" x 1" (3.18 mm x 25.4 mm)
RS-05	PS-05	Solid Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
RS-06	PS-06	Solid Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
RS-07	PS-07	Solid Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
RS-08	PS-08	Solid Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
RS-09	PS-09	Solid Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
RS-10	PS-10	Solid Assorted	Assorted

PACKAGING: Bag of 100

ITEM #	DIMENSIONS
OP-11	0.25" x 0.75" (6.35 mm x 19.05 mm)
OP-12	0.375" x 1" (9.53 mm x 25.4 mm)
OP-13	0.5" x 1" (12.7 mm x 25.4 mm)
OP-14	0.75" x 1" (19.05 mm x 25.4 mm)
OP-15	1" x 1" (25.4 mm x 25.4 mm)
OP-16	1.25" x 1.25" (31.75 mm x 31.75 mm)
OP-17	1.5" x 1.5" (38.1 mm x 38.1 mm)
OP-18	1.75" x 1.25" (44.45 mm x 31.75 mm)
OP-20	Assorted

PACKAGING: Bag of 50


OSTEOTOME PROTECTORS


ITEM #	DIMENSIONS	PACKAGING
RGS-01	1.5" (38 mm)	55 / Bag
RGS-02	2" (51 mm)	45 / Bag
RGS-03	3" (76.2 mm)	40 / Bag

SINGLE-TIP RETRACTOR PROTECTORS



ITEM #	COLOR	SIZE
RGM-01	Pink	Small
RGM-02	Clear	Small
RGM-03	Green	Medium
RGM-04	Clear	Medium

PACKAGING: Bag of 50

WEITLANER RETRACTOR PROTECTORS


ITEM #	COLOR	DIMENSIONS
DG-01	Clear	0.0625" x 0.75" (1.59 mm x 19.05 mm)
DG-02	Blue	0.078" x 0.75" (1.98 mm x 19.05 mm)
DG-03	Green	0.109" x 0.75" (2.78 mm x 19.05 mm)
DG-04	Red	0.125" x 1" (3.18 mm x 25.4 mm)
DG-05	Yellow	0.187" x 1" (4.76 mm x 25.4 mm)
DG-07	Orange	0.062" x 0.375" x 1" (1.59 mm x 9.53 mm x 25.4 mm)
DG-10	Assorted	Assorted

PACKAGING: Bag of 50

PERFORATED DOUBLE PROTECTORS


BATRIK PADLOCKS SECURITY TAGS

Our Padlock Security Tags for instrument containers and general hospital applications are specifically designed for the healthcare market. Our material selection, fabrication process, and engineering ensure consistent quality and high temperature compliance with the various steam sterilization temperatures in worldwide hospital applications.

Available with optional reactive ink indicator dots that provide visual sterility indication. Light blue dots turn brown/dark brown when steam sterilized.



Made of polypropylene



Latex and DEHP Free



Will not melt or soften with any established steam sterilization protocol



Fully open design (45°) facilitates insertion into instrument containers



Easy twist breaking design



Large base for easier handling




Available in different size




Custom Colors and OEM available upon request

AVAILABLE IN THREE FORMATS


700-SERIES



800-SERIES



900-SERIES



PRINTED & CUSTOM NUMBERED PADLOCKS

Available in 700, 800 & 900 SERIES

PRODUCT SELECTION

COLOR	PRINTED WORD
PRINTED PADLOCKS	
WHITE	CLEAN
BLUE	STERILE
RED	REPAIR
YELLOW	SHARPEN
YELLOW	LOANER
BLACK	DISCARD
CUSTOM NUMBERED	
ALL COLORS	up to 8 numbers



CUSTOM NUMBERED
Custom printed padlocks (up to 8 characters)



PRINTED
6 different printed labels available

NOTE: Printed & Custom Numbered Padlocks are available in 700, 800, and 900 Series and can be ordered in 100, 200, and 500 unit quantities. Please contact us for ordering details.

INDICATOR DOTS



REGULAR



700-SERIES PADLOCKS

The 700 series padlocks are perfect for "Aesculap" containers and those with a short, large opening.

PRODUCT SELECTION

REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-701	15-711	15-710
BLUE	15-702	15-712	15-720
GREEN	15-703	15-713	15-730
RED	15-704	15-714	15-740
YELLOW	15-705	15-715	15-750
PINK	15-706	15-716	15-760
ORANGE	15-707	15-717	15-770
BLACK	15-709	15-719	15-790
LIGHT BLUE	15-722	15-732	15-742

WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-701-IND	15-711-IND	15-710-IND
BLUE	15-702-IND	15-712-IND	15-720-IND
GREEN	15-703-IND	15-713-IND	15-730-IND
RED	15-704-IND	15-714-IND	15-740-IND
YELLOW	15-705-IND	15-715-IND	15-750-IND
PINK	15-706-IND	15-716-IND	15-760-IND
ORANGE	15-707-IND	15-717-IND	15-770-IND
BLACK	15-709-IND	15-719-IND	15-790-IND
LIGHT BLUE	15-722-IND	15-732-IND	15-742-IND

Double locking hooks can be inserted in a large guide hole



The grooves on the padlock facilitate handling



INDICATOR DOTS

REGULAR



800-SERIES PADLOCKS

The 800 series padlocks have a longer arm (32 mm) and have a single, thinner hook than the 700 series. They are perfect for almost any type of container, even those with long and thin channels.

PRODUCT SELECTION

REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-801	15-811	15-810
BLUE	15-802	15-812	15-820
GREEN	15-803	15-813	15-830
RED	15-804	15-814	15-840
YELLOW	15-805	15-815	15-850
PINK	15-806	15-816	15-860
ORANGE	15-807	15-817	15-870
BLACK	15-809	15-819	15-890
LIGHT BLUE	15-822	15-832	15-842

WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-801-IND	15-811-IND	15-810-IND
BLUE	15-802-IND	15-812-IND	15-820-IND
GREEN	15-803-IND	15-813-IND	15-830-IND
RED	15-804-IND	15-814-IND	15-840-IND
YELLOW	15-805-IND	15-815-IND	15-850-IND
PINK	15-806-IND	15-816-IND	15-860-IND
ORANGE	15-807-IND	15-817-IND	15-870-IND
BLACK	15-809-IND	15-819-IND	15-890-IND
LIGHT BLUE	15-822-IND	15-832-IND	15-842-IND

Single locking hook can be inserted in a large guide hole




The locking system allows for **easy closure**



INDICATOR DOTS

REGULAR



900-SERIES PADLOCKS

The 900 series padlocks are identical to the 800 series, but with a larger rectangular base which allows for the placement of a 20 mm x 30 mm tracking label.

PRODUCT SELECTION

REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-901	15-911	15-910
BLUE	15-902	15-912	15-920
GREEN	15-903	15-913	15-930
RED	15-904	15-914	15-940
YELLOW	15-905	15-915	15-950
PINK	15-906	15-916	15-960
ORANGE	15-907	15-917	15-970
BLACK	15-909	15-919	15-990
LIGHT BLUE	15-922	15-932	15-942

WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-901-IND	15-911-IND	15-910-IND
BLUE	15-902-IND	15-912-IND	15-920-IND
GREEN	15-903-IND	15-913-IND	15-930-IND
RED	15-904-IND	15-914-IND	15-940-IND
YELLOW	15-905-IND	15-915-IND	15-950-IND
PINK	15-906-IND	15-916-IND	15-960-IND
ORANGE	15-907-IND	15-917-IND	15-970-IND
BLACK	15-909-IND	15-919-IND	15-990-IND
LIGHT BLUE	15-922-IND	15-932-IND	15-942-IND

Single locking hook
for use with smaller diameters



Ideal for adding labels/ indicators indicators or markings to the large rectangular base



GET IN TOUCH

Contact us for bulk orders or
to request your free sample.

USA & Canada Toll Free: 1.800.661.5432

International: +1.514.631.7988

Fax: +1.514.631.9083

Email: info@surgmed.com

