



NEW!



Introducing PandaPatch™: Innovative Adhesive Eye Protection for Safer Patient Care

PandaPatch™ offers an effective alternative to traditional tape for protecting patients' eyes during surgery. These adhesive eye protectors are designed to reduce the risk of corneal injuries and ophthalmic complications for patients under general anesthesia.

With PandaPatch™, you ensure superior patient outcome by preventing avoidable eye damage associated with the use of standard surgical tape.

KEY FEATURES:

- **Ergonomic Design:** Contours comfortably around the eyes.
- **Post-Surgery Friendly:** Particularly effective for children, the elderly, and individuals with sensitive or delicate skin.
- **Latex-Free & Hypoallergenic Adhesive:** Made from medical-grade material, gentle on the skin and designed to minimize irritation.
- **Non-Adhesive Colored Tabs:** Allowing quick and hassle-free application and removal, facilitating effortless monitoring of pupil dilation during surgery.



70-04ADU
Adult Sensitive Eye Protector

HERE'S WHY IT MATTERS:

- During prolonged surgery under general anesthesia, up to 59%¹ of patients may not fully close their eyes, risking corneal abrasions.
- PandaPatch™ gentle adhesive eye protectors ensure eyelids remain safely and comfortably closed, minimizing the risks typically associated with the use of standard surgical tape.



REFERENCES


1. Anson, Jonathan, M.D. "Perioperative Corneal Abrasions: Etiology, Prevention, and Management." Clinical Updates, Penn State University College of Medicine, 2021.




PANDAPATCH™ RANGE:

Three adhesive eye protectors to accommodate different age groups and skin sensitivities. Providing optimal protection for every patient.

- Adult Extra Sensitive Adhesive (Blue Tabs)
- Adult Adhesive (Red Tabs)
- Pediatric Adhesive (Red Tabs)



70-04ADU
Adult Adhesive Eye Protector



70-04PED
Pediatric Adhesive Eye Protector

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
70-02ADU	ADULT EXTRA SENSITIVE ADHESIVE EYE PROTECTORS (BLUE)	2.15" x 1.35"	50 / Box
70-04ADU	ADULT ADHESIVE EYE PROTECTORS (RED)	2.15" x 1.35"	50 / Box
70-04PED	PEDIATRIC ADHESIVE EYE PROTECTORS (RED)	1.85" x 0.93"	50 / Box