

NEW!

TRIFECTA™ BRUSH

INSTRUMENT CLEANING

Introducing the **Trifecta™ Brush**, a versatile **3-in-1 tool** designed to meet all your instrument cleaning needs. This reusable instrument brush is available in two styles: one with **all-nylon** bristles for gentle cleaning and another with nylon bristles on two sides and **stainless steel bristles** on the third, ideal for tackling tougher buildup and stubborn cleaning tasks. This innovative design ensures effective, thorough maintenance for a wide range of instruments.

TRIFECTA™ BRUSH - NYLON

1 NYLON BRISTLES


Gentle, transparent nylon bristles provide an effective yet safe cleaning solution. This side allows you to verify cleanliness easily, without risking damage to delicate surfaces. Ideal for a light, non-invasive clean.

2 PRECISION TIP

Reach those hard-to-clean spots with the precision tip. This small, firm-bristle tip is designed for tight, hard-to-reach areas, such as instrument jaws and box joints. It ensures a deep clean where it's most needed, without the risk of damage or injury to delicate parts.

3 FIRM NYLON BRISTLES

For tougher residues, the firm nylon side delivers the cleaning power you need. With bristles that feel similar to stainless steel but are safe for your instruments, this side easily handles stubborn dirt and buildup.



TRIFECTA™ BRUSH - NYLON & STAINLESS STEEL



→ STAINLESS STEEL BRISTLES ARE MORE DURABLE AND EFFECTIVE FOR CLEANING STUBBORN BUILDUP.

PRODUCT SELECTION

ITEM #	COLOR	DESCRIPTION	LENGTH TOTAL	BRISTLE LENGTH	BRISTLES 1 Ø	BRISTLES 2 Ø	BRISTLES 3 Ø
28-23012	Blue	Triple Nylon Brush	8.66" / 220 mm	1.20" x 1.20" x 0.47" / 30.5 x 30.5 x 12.0 mm	0.47" / 12.0 mm	0.47" / 12.0 mm	0.315" / 8.0 mm
29-23012	Blue	Triple Nylon / Stainless Steel	8.66" / 220 mm	1.20" x 1.20" x 0.47" / 30.5 x 30.5 x 12.0 mm	0.47" / 12.0 mm	0.47" / 12.0 mm	0.47" / 12.0 mm

PACKAGING: Bag of 3 - For a bag of 50, add "-50" following the item code.